



**NEW**

# INTRODUCING THE **SABRE** SCOPE PRE-CLEAN SYSTEM





## WHAT IS SABRE?

Sabre is an all in one channel cleanser and scope sponge designed for the initial flush and wipe of an endoscope, immediately post-procedure, as indicated by most governing bodies. Sabre consists of a two-tiered container, a universal sponge, and a sachet of Dr Weigert Mediclean Forte™ detergent concentrate. All you need to do is add water. The container has a two-tiered level indicating a dilution level of 300ml and 500ml. The dual level indicator allows the Sabre Scope Pre-Clean System to comply with various scope cleaning guidelines.

## HOW DOES SABRE ADHERE TO INTERNATIONAL GUIDELINES?

### World Gastroenterology Organisation

"Pre-cleaning: Clear gloss debris by sucking detergent through the working channel (250ml)"

### European Society of Gastroenterology

"Before the endoscope is detached from the light source and video processor, detergent solution is sucked through the working channel"

### British Society Of Gastroenterology

"Through manual cleaning with a compatible enzymatic detergent, including the brushing and flushing of all accessible endoscope channels, must be undertaken before automatic endoscope disinfection.

This routine must be undertaken during list, between patients and after each patient examination"

FEATURE	BENEFIT
Dr Weigert Mediclean Forte™ detergent	The unique blend of surfactants in Sabre immediately breaks down bio-burden and prevents biofilm from forming
Dedicated and approved for endoscopy use with corrosion inhibitors and anti-foaming agents	Specifically tested for use with endoscopes. The anti-corrosion inhibitors are designed to protect and lubricate scope channel surfaces. Compatibility statements available from Olympus / Pentax / Fujinon for Dr Weigert Mediclean Forte™
Ergonomic handles with grip	Slip free handles to prevent spills and allows for easy transportation
Easy to see how much detergent has been aspirated	Measuring lines between 300ml and 500ml allow the user to monitor aspiration
Easy peel sachet and universal shaped sponge	Detergent and sponge provided in one
Detergent is 100% biodegradable, non-toxic, safe to handle and flush away	Safe for user and environment
Saves on storage space	Packaging design saves on storage, as containers are stackable
No water provided	Using detergent concentrate rather than pre-made detergent is lightweight, reducing transport costs. Less packing material reduces the waste stream

## WHATS INCLUDED?



## ORDER DETAILS

- **Product Code:** SABRE
- **Description:** Sabre Scope Pre-Clean System
- **Box Size:** 100

🏠 **Diagmed Healthcare Ltd**  
Lumley Close, Thirsk Industrial Park, Thirsk  
North Yorkshire, UK, YO7 3TD

☎ Tel: (01845) 526660 📠 Fax: (01845) 522199

✉ Email: [info@diagmed.co.uk](mailto:info@diagmed.co.uk)

🌐 Website: [diagmed.healthcare](http://diagmed.healthcare)



**DIAGMED**  
HEALTHCARE