

# Rules of Conduct

## HaemoCer™ PLUS

Ease of use, DAPI application technique



### DRY

Remove excess pooled blood by suction or gauze – to ensure APH (Absorbable Polysaccharide Haemostat) particles are activated at the wound site



### WET

The APH particles will be activated on the top of the pooled blood only and will not assist clot formation at the wound site

### APPLY

Apply HaemoCer™ PLUS immediately and extensively to the wound site ensuring the entire wound perimeter is fully and swiftly covered



### APPLY

Avoid applicator contact with blood

Don't construct small heaps or slowly and gradually apply APH particles

### PRESS

Apply gentle compression over the wound using **DRY** gauze without any interaction



### PRESS

No twisting or turning with compression gauze

Using a wet gauze will activate APH particles from the top and risks pulling away the gelled matrix


### IRRIGATE

Irrigate thoroughly prior to gently removing the gauze to ensure the gelled matrix remains at the wound site and is not torn away with the gauze



### TEAR OFF

Don't disturb the gelled matrix after haemostasis has been achieved



made in  
Germany

# HaemoCer<sup>TM</sup> PLUS

## Haemostatic Powder

HaemoCer<sup>TM</sup> PLUS – Absorbable Polysaccharide Haemostat (APH) engineering incorporates a sophisticated, plant-based polymer crosslinking (Polysaccharide Ultra-hydrophilic Resorbable Engineering) PURE<sup>TM</sup> processing that creates ultra-hydrophilic, biocompatible particles.

Utilizing a dual mode of action upon contact with blood HaemoCer<sup>TM</sup> PLUS firstly enhances the natural clotting cascade by rapidly dehydrating the blood and accelerating the concentration of platelets, red blood cells and coagulation proteins at the bleeding site.

The second mode of action of HaemoCer<sup>TM</sup> PLUS is the formation of a robust gelled matrix that adheres to and forms a mechanical barrier to further bleeding.

HaemoCer<sup>TM</sup> PLUS contains no animal or human components and is rapidly resorbed by amylase within 48 hours without leaving any residues behind.



100% resorbable



resorbed within 48h



Plant based



Ease of use

# BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Str. 36c · 95447 Bayreuth · Germany · Tel. +49 (0)921 78 77 70 0

Fax. +49 (0)921 78 77 70 79 · info@biocer-gmbh.de · www.biocer-gmbh.de

