Roth Net® retriever – 360

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Roth Net®-Bergungsinstrument – 360

Retino di recupero Roth Net® – 360

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*listening... and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

Not made with natural rubber latex.

Intended Use:

The Roth Net® retriever – 360 product line is intended to be used to retrieve excised polyps, tissue samples, foreign bodies and calculi during flexible and rigid endoscopy procedures.

Description	Product Number	Sheath Diameter	Minimum Scope Channel required	Sheath Length	Rotation	Units/Box
Roth Net® retriever – 360 (sterile)	00711198	2.5mm	2.8mm	230cm	Yes	10

Warnings and Precautions:

- The endoscopic retrieval of foreign objects, food bolus, or polyps should only be performed by persons having adequate training and familiarity with endoscopic techniques.
- Consult the medical literature relative to techniques, complications and hazards prior to the performance of any endoscopic procedure.
- If reference information for endoscopic techniques utilizing the Roth Net® retriever product line is desired/needed, contact US Endoscopy or refer to www.usendoscopy.com. Case reports and videos describing and demonstrating techniques are available.
- Care should be exercised when grasping the targeted object, material, or polyp to be retrieved to avoid inadvertently grasping tissue or organs not intended for retrieval.
- **Keep gentle traction on the device during retrieval so that the retrieved object does not become loosened, dislodged, or aspirated into the trachea.**
- **If resistance to insertion is encountered, reduce the angulation (or lower the forceps elevator if applicable) until the instrument passes smoothly.**
- **Short strokes, 1"-1.5" (2.5cm – 3.8cm) in length, are recommended throughout device passage to avoid sheath kinking.**
- The technique for using the Roth Net® retriever – 360, as described in any document or publication, is at the discretion of the medical personnel attempting the retrieval.
- This device is not recommended for the retrieval of sharp foreign objects.
- This device is not recommended for use in side viewing endoscopes.
- The following conditions may cause the device to function improperly:
 - **Advancing the handle to the open position with too much speed or force.**
 - **Attempting to pass or open the device in an extremely articulated endoscope.**
 - **Actuating the device in an extremely coiled position.**
 - **Actuating the device when the handle is at an acute angle in relation to the sheath.**
- Do not attempt to reuse, reprocess, refurbish, remanufacture or resterilize this device. US Endoscopy did not design this device nor is it intended to be reused, reprocessed, refurbished, remanufactured, or resterilized. Performing such activities on this disposable medical device presents a safety risk to patients (i.e. compromised device integrity, cross-contamination, infection).
- This device has not been designed nor is it intended to be used with electrosurgical devices or their accessories.

Contraindications:

- Contraindications include, but are not limited to, those specific to any endoscopic procedure.

Prior to Use:

1. Read the entire *Instructions for Use*.
2. Ensure the accessory channel diameter of the endoscope is compatible with the outer diameter of the device catheter prior to insertion (**refer to chart above**).
3. Remove the device from the package.
4. Inspect and familiarize yourself with the device and review the diagrams. If there is evidence of damage to the device (i.e. cracked/kinked catheter, deformed net), **do not use this product and contact your local US Endoscopy representative**.
5. Uncoil the entire device and drape in a "U" shaped configuration.
6. Hold the Proximal End of the catheter in one hand and the Distal End in the opposite hand.
7. Close (retract) and open (deploy) the handle several times and observe that the device functions properly.
8. Close (retract) the handle until the net is completely repacked into its catheter (See FIG. 1).

Directions for Use:

1. Once the object has been endoscopically identified, advance the retracted/closed device into the accessory channel of the endoscope using short strokes (1" – 1.5" in length) until the distal end of the sheath is endoscopically visualized.
2. Advance the distal end of the sheath slightly past the object, bolus, or polyp to be retrieved. **Avoid blindly passing the device past any object, bolus, or polyp if the entire lumen is blocked.**
3. Open (deploy) the device by advancing the handle forward until it stops. Confirm the net is fully open via endoscopic observation (see FIG. 2).
4. The following conditions may cause the device to function improperly:
 - Advancing the handle to the open position with too much speed or force.
 - Attempting to pass or open the device in an extremely articulated endoscope.
 - Actuating the device in an extremely coiled position.
 - Actuating the device when the handle is at an acute angle in relation to the sheath.

5. Endoscopically position the Roth Net® retriever – 360 over the object, bolus, or polyp.

Note: For best results, the object should be centered in the net prior to net closure, however food bolus retrieval techniques may vary.

6. Improve positioning of the net by rotating the handle in the desired direction. (See FIG 3).

Note: Rotation is best achieved with the net open and centered in the lumen, not lying against the mucosal wall.

7. Close (retract) the handle until resistance is felt and the object, bolus, or polyp is entrapped in the net (see FIG. 4).

8. Keep the device closed by applying continuous/steady traction to the handle.

Note: Excessive closing force may damage the device's wire or cause the net to rupture or tear.

Exerting excessive "bending" pressure on the open device may disfigure the wire shape.

9. Grasp and slowly pull the proximal end of the device to bring the collected object close to the tip of the endoscope making sure not to obscure the endoscopic view.

10. Endoscopic observation is necessary during removal as to not lose sight of the object, bolus, polyp, or surrounding tissue during removal.

11. When the object is secure in the net, withdraw the object together with the endoscope from the patient.

12. Do not retract the device and the object into the endoscope, this could cause damage to the endoscope and the device.

13. Once the endoscope and the device have been withdrawn, the retrieved object, bolus, or polyp can be removed from the net by advancing the handle forward to open the device (see FIG. 2).

Note: Rinse net in sterile water to release retrieved object. Do not use fingers to remove object. This could damage the net.

14. Specimens should be prepared for histologic evaluation per institutional guidelines.

15. Confirm that the device continues to open and close smoothly **prior to reinsertion into the patient** (if removing additional objects is required).

Note: Rinsing the net in sterile water may be required if debris remains in net.

16. Close (retract) the handle until the net is completely repacked into its sheath (see FIG. 1).

17. Repeat steps 1-16 as needed.

18. Once the object or objects have been retrieved and removed from net, close (retract) the net into the catheter (see FIG. 1) and remove the device from the endoscope.

FIG. 1


FIG. 2


FIG. 3



FIG. 4

Product Disposal:

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Issued Date: January 2019

Warning:

An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact US Endoscopy to determine if additional information is available.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

US Endoscopy, a wholly owned subsidiary of STERIS Corporation.

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Listening...and delivering solutions®, le logo US Endoscopy listening...and delivering solutions, ainsi que toutes les marques portant le signe ® ou ™ sont des marques déposées auprès de l'U.S. Patent and Trademark Office ou sont des marques commerciales détenues par US Endoscopy.
Listening...and delivering solutions®, das „listening...and delivering solutions“-Design und alle mit ® oder ™ versehenen Marken sind beim United States Patent and Trademark Office (Patentamt der Vereinigten Staaten) registriert oder befinden sich im Besitz von US Endoscopy.
Listening...and delivering solutions®, o design US Endoscopy listening...and delivering solutions e todas as marcas assinaladas com ® ou ™ estão registadas no Departamento de Registos de Patentes dos EUA (U.S. Patent and Trademark Office) ou são marcas comerciais detidas pela US Endoscopy.
Listening...and delivering solutions®, designet af US Endoscopy listening...and delivering solutions, og alle mærker markeret med ® eller ™ er registreret hos U.S. Patent and Trademark Office eller varemærker ejet af US Endoscopy.
Listening...and delivering solutions[®], het ontwerp Listening...and delivering solutions van US Endoscopy en alle aanduidingen met [®] of [™] zijn gedeponeerd bij het Patent and Trademark Office in de Verenigde Staten of zijn handelsmerken van US Endoscopy.
Listening...and delivering solutions®, US Endoscopy listening...and delivering solutions 디자인, 그리고 ® 또는 ™으로 표시된 모든 상표는 미국 특허청(U.S. Patent and Trademark Office)에 등록되어 있거나 US Endoscopy가 소유하고 있는 상표입니다.