

Carr-Locke injection needle



listening...and delivering solutions®

Trusted performance...

World-renowned for its reliability and consistency, the **Carr-Locke injection needle** is designed to operate in the most tortuous conditions.

- **Stainless steel, Teflon coated, spring sheath** to virtually eliminate sheath kinking, while offering exceptional functionality and needle performance during challenging procedures
- **Ultra smooth distal metal hub** aids in tamponade and is gentle on friable tissue, helping to prevent sheath perforations
- **Luer-lock, spring-loaded handle** automatically retracts the needle, reducing the risk of inadvertent needle sticks or scope damage

Consistently deploys in retroflexed position

| product number | sheath diameter (mm) | length (cm) | needle projection (mm) | needle gauge | units/box |
|---------------------|----------------------|-------------|------------------------|--------------|-----------|
| BX00711811 | 2.5 | 230 | 5 | 25 | 5 |
| BX00711812 | 2.5 | 230 | 5 | 23 | 5 |
| BX00711813 (Pentax) | 2.5 | 230 | 5 | 25 | 5 |
| BX00711814 | 2.5 | 230 | 4 | 23 | 5 |
| BX00711822 | 1.8 | 230 | 5 | 25 | 5 |
| BX00711823 | 2.5 | 230 | 5 | 25 | 10 |
| BX00711824 | 2.5 | 230 | 5 | 23 | 10 |

Distributed by:

US endoscopy

*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

760964 Rev. B

US Endoscopy is a subsidiary of STERIS Corporation.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

PENTAX® is a registered trademark of Hoya Corporation

Teflon® is a registered trademark of E.I. duPont de Nemours and Company as its affiliates