



## A variety of needles for a variety of needs...

US Endoscopy features a full line of injection needles to suit nearly all of your injectio varices, ulcers, and tattooing.

## Carr-Locke injection needle

Designed to consistently operate in the most tortuous conditions

now available in 23 gauge

Carr-Locke injection needle

## features and benefits chart

feature	benefit	Carr-Locke needle	Articulator™ needle	Vari-Safe <sup>™</sup> needle
Teflon coated spring sheath and pre-load with stop at distal end	<ul> <li>Ensures consistent and full length needle projection every time</li> <li>Allows for smooth insertion down endoscope channel</li> </ul>	<b>✓</b>		
Stainless steel spring sheath	Eliminates sheath kinking especially during challenging procedures	<b>✓</b>	<b>✓</b>	
Ultra smooth distal metal hub	<ul><li>Controls depth of needle insertion; prevents needle perforations</li><li>Aids in tamponade</li></ul>	<b>√</b>	<b>✓</b>	
Luer lock, spring loaded handle	Automatically retracts needle to reduce risk of inadvertent needle sticks and/or scope damage	<b>√</b>	<b>✓</b>	<b>✓</b>
Requires less than 1cc of fluid to prime (most needles require 1.5-2.0cc)	<ul> <li>Less expensive, particularly when using costly injection agents</li> </ul>	<b>✓</b>	<b>✓</b>	<b>✓</b>

n needs - from saline assisted polypectomy and enteroscopy to esophageal

Articulator™ injection needle

Ideal for everyday use and dependable for your most immediate needs

> Articulator™ injection needle

Vari-Safe<sup>™</sup> injection needle Cost effective and reliable

Vari-Safe™ injection needle

For more information on the injection needles, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	needle projection (mm)	needle gauge	units/box
BX00711803	Articulator™ injection needle	2.5	160	4	25	5
BX00711804	Articulator™ injection needle	2.5	160	5	25	5
BX00711807	Articulator™ injection needle	2.5	230	5	25	5
BX00711808*	Articulator™ injection needle	2.5	350	5	25	5
BX00711810	Articulator™ injection needle	2.5	230	5	25	10
BX00711811	Carr-Locke injection needle	2.5	230	5	25	5
BX00711812	Carr-Locke injection needle	2.5	230	5	23	5
BX00711813 (Pentax compatible)	Carr-Locke injection needle	2.5	230	5	25	5
BX00711814	Carr-Locke injection needle	2.5	230	4	23	5
BX00711819	Vari-Safe™ injection needle	2.5	230	5	23	1
BX00711822	Carr-Locke injection needle	1.8	230	5	25	5
BX00711823	Carr-Locke injection needle	2.5	230	5	25	10
BX00711824	Carr-Locke injection needle	2.5	230	5	23	10

Distributed by:



5976 Heisley Road Mentor, OH 44060 USA phone +1 440 / 639-4494 fax +1 440 / 639-4495 email global@usendoscopy.com

www.usendoscopy.com