

trusted security

**Trusted
Security**



us endoscopy

listening...and delivering solutions[®]

setting the standard

Trusted Security

US Endoscopy and *Bennett Roth, M.D.* partnered to develop a device that revolutionized the way the GI community approaches polyp and foreign body retrieval.

There are many styles of Roth Net[®] retrievers to support a variety of retrieval needs. Securing large and small foreign bodies, preserving tissue specimens - the Roth Net[®] retriever is a device you can count on for consistent, superior performance.

*Roth Net[®] retriever – a device you trust...
from a company you can depend on*



foreign body

foreign body retrieval

The **Roth Net®** and **Roth Net® Platinum®** retrievers are ideal for...

*removing round, blunt or otherwise hard to grasp/
hard to hold foreign bodies*

*securely extracting foreign bodies past the
cricopharyngeus*



Bottle cap retrieval
with Roth Net®
retriever

*Engineered with a strong,
reinforced net for extra
security*



Food bolus removal
with Roth Net®
Platinum® retriever

*Unique octagonal design to
maintain shape and open to
capacity across multiple passes*

confidence

tissue retrieval

tissue retrieval

The **Roth Net® retriever** and **Roth Net® Platinum® retriever** are specifically designed to...

retrieve polyps and polyp fragments quickly and easily without losing them

preserve the morphology of specimens for histopathologic review

allow continued examination with easy visibility through a more transparent net material



Roth Net® Platinum® retrieval of large polyp



Roth Net® retrieval of tissue sample



Roth Net® retriever

trusted

NEW!


Roth Net® retriever - mini polyp

Featuring a smaller net design, offering enhanced maneuverability and visualization during retrieval within the small confines of the colon.

round snare shape improves maneuverability and allows for a clear view during retrieval

contains memory points in the snare that enable net to open to full capacity

maintains a consistent shape during multiple uses



*Roth Net® retriever -
mini polyp actual size
(2.5 x 4.5cm)*



time savings

www.usendoscopy.com

specialty retrieval

specialty retrieval

Roth Net[®] retrievers are specifically designed to address...

unique retrieval needs including small, medium, and large foreign bodies

- **mini** retriever for smaller foreign bodies and smaller channel endoscopes
- **maxi** retriever for larger foreign body removal (AA battery)
- **enteroscopic** retriever for excised polyp retrieval in the small bowel



Multiple coin removal with Roth Net[®] retriever



AA battery retrieval with Roth Net[®] retriever - maxi

mini retriever - now available sterile

peace of mind

product information

part number	description	sheath diameter (mm)	length (cm)	net size (cm)	qty/box
foreign body retrieval					
BX00711050	Roth Net® retriever - standard	2.5	230	3 x 6	5
BX00711052	Roth Net® retriever - standard	2.5	230	3 x 6	10
BX00711053	Roth Net® retriever - standard	2.5	160	3 x 6	5
BX00711155	Roth Net® Platinum® retriever - food bolus	2.5	160	4 x 5.5	5
BX00711162*	Roth Net® retriever - standard (sterile)	2.5	230	3 x 6	10
BX00715050	Roth Net® Platinum® retriever - universal	2.5	230	4 x 5.5	5
BX00715060*	Roth Net® Platinum® retriever- universal (sterile)	2.5	230	4 x 5.5	5
tissue retrieval					
BX00711150	Roth Net® retriever - polyp	2.5	230	3 x 6	5
BX00711151	Roth Net® retriever - polyp	2.5	230	3 x 6	10
BX00715150	Roth Net® Platinum® retriever - polyp	3.0	230	4 x 5.5	5
BX00715160*	Roth Net® Platinum® retriever - polyp (sterile)	3.0	230	4 x 5.5	5
specialty retrieval					
BX00711057	Roth Net® retriever - mini	1.8	160	2 x 4.5	5
BX00711059	Roth Net® retriever - maxi	3.0	160	4 x 8	5
BX00711152	Roth Net® retriever - enteroscopic	2.5	350	3 x 6	5
NEW! BX00711158	Roth Net® retriever - mini polyp	2.8	230	2.5 x 4.5	10
NEW! BX00711159	Roth Net® retriever - mini polyp (sterile)	2.8	230	2.5 x 4.5	5
BX00711165	Roth Net® retriever - maxi (sterile)	3.0	230	4 x 8	5
NEW! BX00711168	Roth Net® retriever - mini (sterile)	1.8	160	2 x 4.5	5

*not currently available for sale in the United States

For customers outside of the United States, please contact US Endoscopy corporate headquarters at global@usendoscopy.com for distributor information.



listening... and delivering solutions™



Roth Net® and Roth Net® Platinum® is a registered trademark of US Endoscopy, U.S. Patent Numbers 6,814,739 and 5,643,283, patents pending.
US Endoscopy is a subsidiary of STERIS Corporation.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and, unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226

Follow us:    

www.usendoscopy.com