



Traxtion®
snare



*uniquely designed,
stiff .018" diameter
braided wire*

A stiff braided snare that grips the mucosa...

The unique braided wire of the Traxtion® snare provides exceptional grip of the mucosa and is ideal for en bloc resection of large flat lesions.



memory points for better grip of the mucosa

The Traxtion® hexagonal snare features...

- A uniquely designed, braided snare wire that makes it stiffer than other hexagonal snares
- Memory points for better grip of the mucosa
- The ability to retain its shape and integrity with multiple resections

For more information on the Traxtion® snare, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	snare size (width x length)	active cord style	unit per box
BX00711113	Traxtion® snare	2.4	230	25mm x 40mm	Olympus	10

Distributed by:

us endoscopy

listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

760375 Rev. C

US Endoscopy is a subsidiary of STERIS Corporation.
© 2015 US Endoscopy or its affiliates. All rights reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

Follow us:    