

eTrap[®] polyp trap



The ideal polyp trap designed with you in mind...

Inspired by nurses, the eTrap[®] polyp trap allows for quick and accurate retrieval of multiple specimens, while safeguarding clinicians from unnecessary exposure to biomaterial.

The eTrap[®] polyp trap features...

- two removable specimen strainers designed for retrieval of multiple polyps with uninterrupted suction
- clear magnifying window for easy visualization of collected specimen
- simple transport and preservation of specimens



specimen captured in eTrap[®] polyp trap

For more information on the eTrap[®] polyp trap, contact your local representative today.

product number	description	tubing	units/box
BX00711099	polyp trap	~ 15 cm long	25 units/box 2 strainers/unit

Distributed by:



listening...and delivering solutions[®]

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

US Endoscopy is a subsidiary of STERIS Corporation.
© US Endoscopy or its affiliates. All rights reserved.
Listening...and delivering solutions[®], the US Endoscopy listening...and delivering solutions design, and all marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

Follow us:    