

DuoSwift™ combination squeegee brush



Combining the cleaning benefits of a bristled brush and squeegee into one device...

The DuoSwift™ combination squeegee brush features a tapered nylon brush on the leading end and a series of squeegees on the trailing end to effectively remove debris from flexible endoscopes while adhering to industry guidelines.



The DuoSwift™ combination squeegee brush features a...

- bristled brush to help loosen and remove bioburden with friction
- squeegee to pull liquid and remaining residue through the endoscope channel
- distinct design to help maximize cleaning touches within compatible endoscope channels

For more information on the DuoSwift™ combination squeegee brush, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	bristle diameter (mm)	squeegee diameter (mm)	channel size (mm)	units/box
EX00711618	DuoSwift™ combination squeegee brush	1.7	230	5-7	2.5-5.2	2.8-7.0	50
EX00711619	DuoSwift™ combination squeegee brush with valve/control head cleaning brush	1.7	230	5-7 10 (valve/ control head)	2.5-5.2 (squeegee)	2.8-7.0	50

Distributed by:



listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

US Endoscopy is a subsidiary of STERIS Corporation.

©2016 US Endoscopy. All Rights Reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

Follow us:    

761177 Rev. A