

QuickCatch™
in-line suction polyp trap

Suction Polyp Trap™
four chamber trap

**SUCTION
POLYP TRAP™**



In-line polyp retrieval solutions...

The Quick Catch™ in-line suction polyp trap and Suction Polyp Trap™ four chamber trap allow for easy in-line installation and feature clear plastic construction allowing for quick visual verification of polyp specimens.



The **Quick Catch™ in-line suction polyp trap** is a lower-cost option featuring a...

- Unique design - offering a chamber with flat surfaces allowing for easy gripping and handling
- Chamber size that easily fits into most specimen jars, reducing polyp handling



The **Suction Polyp Trap™ four chamber trap** is a safe and effective option...

- Featuring four capture chambers allowing for accurate identification and localization of multiple specimens
- Designed to filter specimens helping to minimize time and exposure to bodily fluids

For more information on these products, contact US Endoscopy today.

product number	description	units/box
710201	Quick Catch™ in-line suction polyp trap	50
710202	Suction Polyp Trap™ four chamber trap	12

Distributed by:



*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

US Endoscopy is a subsidiary of STERIS Corporation.

© 2016 US Endoscopy or its affiliates. All rights reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

Follow us:    