




pH Neutral Scope Pre-Clean *Now Available*

MAKING A DIFFERENCE TO HEALTH

SABRE pH is a safe, reliable, compatible and neutral detergent specifically designed for the bedside preclean of flexible endoscopes. Just add water to dilute the concentrated detergent to the preferred volume as outlined in scope manufacturers guidelines.



FEATURES

- Product contains biodegradable organic materials and is gentle on internal scope channels
- Compatible with all metal, rubber, plastic, glass and surgical instruments
- Contains stabilising enzyme which remains inactive during the bedside preclean process

HOW DOES SABRE PH IMPROVE PERFORMANCE?

Its rapid contact time creates a biological catalyst to increase the chemical activity within the detergent solution. This reaction helps to break the bioburden into smaller molecules for easier removal within the decontamination process.

¹"Endoscope channels should be flushed immediately following the procedure and the external surfaces socially cleaned" ¹National Endoscopy Programme - Decontamination Standards for Flexible Endoscopes, March 2008

²"Pre-cleaning.... Flush all channels with a low foaming neutral detergent until runs clear"

²BSG Guidance on Decontamination of Equipment for Gastrointestinal Endoscopy: 2014 Edition

ORDERING INFORMATION

PRODUCT CODE	DESCRIPTION	BOX SIZE
SABRE PH	PH NEUTRAL SCOPE PRE-CLEAN SYSTEM	100

🏠 **Diagmed Healthcare Ltd**
Lumley Close, Thirsk Industrial Park, Thirsk
North Yorkshire, UK, YO7 3TD

☎ Tel: (01845) 526660 📠 Fax: (01845) 522199

✉ Email: info@diagmed.co.uk

🌐 Website: diagmed.healthcare



DIAGMED
HEALTHCARE