

Talon™ grasping device



US endoscopy

listening... and delivering solutions®

Offering an extra hand in food impaction retrieval...

The **Talon™ grasping device** is specifically designed to provide the firm hold required for challenging food bolus and foreign body retrieval cases.

Featuring **4 stainless steel prongs** with **inward-facing atraumatic hooks**, the Talon™ grasping device was designed to mimic the grasping capabilities of a sea hawk.

Uniquely engineered, the device **helps grab into the meat impaction** or around the foreign body providing the strength required for a swift retrieval.

The Talon™ grasping device is designed to remove:

- *meat impactions*
- *food in the esophagus due to achalasia*
- *foreign bodies*



inward-facing atraumatic hooks

For more information on the Talon™ grasping device, contact US Endoscopy today.

product number	description	sheath diameter	sheath length	units/box
00711175	Talon™ grasping device	2.5mm	160cm	5

Distributed by:



*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email global@usendoscopy.com

www.usendoscopy.com

760352 Rev. C