

# bite blocks



**us** endoscopy

*listening... and delivering solutions®*

## ***A broad offering to fit all patient needs...***

Our extensive line of bite blocks offers a variety of options in sizes, designs, and models to ensure the best for the patient and type of procedure.

### **Standard**

The sturdy, yet forgiving material of the standard bite blocks protects delicate instruments and provides comfort for patients.



### **maxi bite block / 60FR**

- available with a VELCRO® Brand fastener, latex-free or woven strap or without a strap
- with dental retention rim



### **standard bite block / 48FR**

- available with a VELCRO® Brand fastener, latex-free or woven strap, or without a strap
- with dental retention rim



### **intermediate bite block / 48FR**

- latex-free strap
- without dental retention rim

## Specialty

From a pediatric bite block to a purpose-built O<sub>2</sub> delivery bite block, our specialty bite blocks are designed to fit the unique needs of patients.



### E-Z-Guard™ bite blocks

- sturdy, inner-core construction – protecting your device or instrument
- soft, foam exterior – providing greater comfort for patients
- individually wrapped package



### pediatric bite block / 38FR

- latex-free strap
- without dental retention rim



### Vasto™ bite block / 72FR

- featuring larger side ports for scope repositioning and suction access
- latex-free strap
- with dental retention rim



### Respa™ O<sub>2</sub> delivery bite block / 60FR

- can be used to house the nasal cannula during procedures to easily deliver O<sub>2</sub> to the patient orally
- latex-free strap
- without dental retention rim

For more information, contact US Endoscopy today.

product number	description	opening size	dental retention rim	units/box
<b>Standard</b>				
BX00712801	standard - latex-free strap	16mm / 48 FR	yes	100
BX00712802	standard - no strap	16mm / 48 FR	yes	100
BX00712803	standard - VELCRO® Brand fastener strap	16mm / 48 FR	yes	100
BX00712804	maxi - latex-free strap	20mm / 60 FR	yes	100
BX00712805	maxi - no strap	20mm / 60 FR	yes	100
BX00712806	maxi - VELCRO® Brand fastener strap	20mm / 60 FR	yes	100
BX00712809	intermediate - latex-free strap	16mm / 48 FR	no	100
BX00712812	maxi - woven strap	20mm / 60 FR	yes	100
<b>Specialty</b>				
BX00712808	pediatric - latex-free strap	12.6mm / 38 FR	no	25
BX00712810	Vasto™ bite block - latex-free strap	24mm / 72 FR	yes	100
BX00712811	Respa™ O <sub>2</sub> delivery bite block	20mm / 60 FR	yes	100
710101	E-Z-Guard™ bite block with foam	48 FR	-	24
710103	E-Z-Guard™ bite block with foam and strap	48 FR	-	24
710106	E-Z-Guard™ pediatric bite block with foam and strap	36 FR	-	24

All E-Z-Guard™ bite blocks are manufactured not to include latex.



5976 Heisley Road  
 Mentor, OH 44060  
 phone (1) 440 / 639-4494  
 fax (1) 440 / 639-4495  
 customer service (1) 800 / 769-8226

[www.usendoscopy.com](http://www.usendoscopy.com)

VELCRO® is a registered trademark of Velcro BVBA. Used with permission.  
 US Endoscopy is a subsidiary of STERIS Corporation.

© 2017 US Endoscopy or its affiliates. All rights reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

Follow us:    