

polyp trap retrieval solutions

listening...and delivering solutions[®]

US Endoscopy is proud to offer a complete portfolio of ***polyp trap retrieval solutions*** that address a variety of Customer preferences and clinical needs.


From the innovative ***eTrap[®] polyp trap*** designed with nurses in mind, to the economical ***Baleen[™] polyp trap***, our portfolio includes a variety of designs to accommodate you and your staff.

eTrap[®] polyp trap

Inspired by nurses, the **eTrap[®] polyp trap** is our most popular solution allowing for quick and accurate retrieval of multiple specimens, while safeguarding clinicians from unnecessary exposure to biomaterial.

The eTrap[®] polyp trap...

- *includes 2 removable strainer trays*
- *allows for collection of multiple polyps with uninterrupted suction*
- *features a clear magnifying window for easy visualization of collected polyps*
- *offers simple transport and preservation of specimens*



specimen captured in eTrap[®] polyp trap



Suction Polyp Trap™ four chamber trap

The **Suction Polyp Trap™ four chamber trap** helps reduce procedure time by eliminating the need for hand-screening and filtering, allowing for quick retrieval of polyps.

The Suction Polyp Trap™ four chamber trap...

- *features four capture chambers allowing for accurate identification and localization of multiple specimens*
- *offers a rotatable lid for efficient polyp retrieval*
- *is a safe and effective option, designed to filter specimens, helping to minimize time and exposure to bodily fluids*



**SUCTION
POLYP TRAP™**



Quick Catch™ in-line suction polyp trap

The **Quick Catch™ in-line suction polyp trap** allows for easy in-line installation and features clear plastic construction allowing for quick visual verification of polyp specimens.


Quick Catch™
In-Line Suction
Polyp Trap

The Quick Catch™ in-line suction polyp trap...

- offers a chamber with flat surfaces allowing for easy gripping and handling
- features a luer connector to flush polyps from the screen with a syringe



luer lock connector



Baleen™ polyp trap

The **Baleen™ polyp trap** reduces exposure to biomaterial and is an economical alternative to other polyp trap solutions.



The Baleen™ polyp trap...

- *easily attaches to suction tubing to act as an in-line filter collecting the specimen*
- *offers a sturdy, durable netting material that helps to resist tearing*



Baleen™ polyp trap inserted into suction tubing.



Baleen™ polyp trap in place and ready for use.



Specimen collected in netting of Baleen™ polyp trap.

product information

product number	description	units/box
BX00711099	eTrap® polyp trap	25
710202	Suction Polyp Trap™ four chamber trap	12
710201	Quick Catch™ in-line suction polyp trap	50
NBPT520	Baleen™ polyp trap	100

For more information on our complete offering of polyp trap solutions, visit www.usendoscopy.com/polypectomy today.



Follow us:



US Endoscopy is a subsidiary of STERIS Corporation.

© 2018 US Endoscopy or its affiliates. All rights reserved.

All marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.

5976 Heisley Road
Mentor, OH 44060

phone (1) 440 / 639-4494

fax (1) 440 / 639-4495

customer service (1) 800 / 769-8226

www.usendoscopy.com

761490 Rev. B