



DIAGMED
HEALTHCARE



NEW

INTRODUCING  **AQUAVIPER**
THE CLEARER WAY TO HAEMOSTASIS

THE LATEST ADDITION TO THE VIPER RANGE

The AquaViper Haemoclclip is a unique flushing catheter to assist with clear visualisation when clipping, with one-to-one rotation for precise placement, even in the most tortuous positions



MAKING A DIFFERENCE TO HEALTH



AQUAVIPER

HAEMOCLIP

FEATURE	BENEFIT
2-way 360° Jaw Rotation	One-to-one rotation to ensure precise positioning and placement
Integral Flushing Port	Can be used in conjunction with external water pump
Unlimited Repositioning	Avoids wasted clips and shortens procedure time
135° Opening Angle	Allows for greater flexibility of placement in difficult areas of the GI Tract
Short Clip Length	Easier to place multiple clips without overlap. Easier deployment and clip detachment - One Click and It's Gone!
Compatibility with a wide range of scopes	230cm catheter length. Flexible catheter allows placement in retroflex positions and via a side-viewing scope
MRI Compatible	Nitinol clip ensures MRI compatibility (Conditional to 3 TESLA)



ORDERING INFORMATION

PRODUCT CODE	DESCRIPTION	CLIP OPENING (MM)	BOX SIZE
AQUA VIPER	AQUAVIPER Haemoclip	11	5

Diagmed Healthcare Ltd

Lumley Close, Thirsk Industrial Park, Thirsk
North Yorkshire, UK, YO7 3TD

Tel: (01845) 526660 Fax: (01845) 522199

Email: info@diagmed.co.uk

Website: diagmed.healthcare



DIAGMED

HEALTHCARE