

**Guardus<sup>®</sup>**  
**overtube**



# Offering protection and peace of mind...


Specifically developed for enhanced protection and safety, the Guardus® overtube system features an exclusive tapered tip design to protect against mucosal damage, while its unique insufflation cap helps maintain insufflation and guard against leakage of body fluids.

**The Guardus® overtube is designed to assist with the following endoscopic needs:**

- foreign body/sharp object extraction
  - stent removal
  - upper endoscopic mucosal resection
- minimizing the risk of aspiration
  - gastroparesis and achalasia
  - active upper GI bleed
- procedures requiring multiple scope insertions
  - food bolus procedures



insufflation cap



inner tube with tapered tip

overtube inner diameter measures 16.7mm

For more information on the Guardus® overtube, contact US Endoscopy today.

product number	description	scope OD (mm)	length (cm)	overtube OD (mm)	overtube ID (mm)	units/box
00711146	esophageal length/standard	8.6-10.0	25	19.5	16.7	3
00711147	esophageal length/therapeutic	10.0-11.7	25	19.5	16.7	3
00711148	gastric length/standard	8.6-10.0	50	19.5	16.7	3
00711149	gastric length/therapeutic	10.0-11.7	50	19.5	16.7	3

Distributed by:



listening...and delivering solutions®

5976 Heisley Road  
Mentor, OH 44060 USA

phone +1 440 / 639,4494

fax +1 440 / 639,4495

email [global@usendoscopy.com](mailto:global@usendoscopy.com)

[www.usendoscopy.com](http://www.usendoscopy.com)

760339 Rev. A