

# Raptor® grasping device



# Offering the precision and strength needed during foreign body retrieval...



## *The Raptor® device:*

- combines alligator and rat tooth capabilities into one device, increasing the gripping ability for a stronger hold
- eliminates reprocessing and ensures the device performs consistently each time it is used
- operates in a retroflexed position and in side-viewing endoscopes

## *NEW Raptor® grasping device - mini:*

- features a smaller hybrid jaw and catheter option for tight or torqued positions
- is specially designed for procedures requiring an ultra-slim endoscope
- offers the firm grip needed to retrieve smaller foreign bodies and stents



For more information on the Raptor® grasping device, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	sterile	units/box
00711177	Raptor® grasping device	2.4	230	yes	5
00711178	Raptor® grasping device - mini	1.9	200	yes	5

Distributed by:

**US endoscopy**

*listening...and delivering solutions\**

5976 Heisley Road  
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email [global@usendoscopy.com](mailto:global@usendoscopy.com)

[www.usendoscopy.com](http://www.usendoscopy.com)

760750 Rev. E