

iSnare[®] system



“The advantage of the combined snare over the braided snare, as seen in our study, is likely multifactorial. This study provides evidence that the integrated needle-snare may be superior to the snare alone for the removal of large, flat polyps*.”

**T. Woodward, J. Crook, M. Raimondo, M. Wallace. "Improving complete EMR of colorectal neoplasia: a randomized trial comparing snares and injectate in the resection of large sessile colon polyps." GIE Journal online, 2015.*

Combining a needle and snare into one device to simplify advanced polyp resection...

The iSnare® system is a family of unique, multi-functional devices ideal for minimizing procedure time and optimizing clinical outcomes during submucosal injection polypectomy procedures including EMR (endoscopic mucosal resection) and SAP (saline assisted polypectomy).

The iSnare® system

- Features alternating injection needle and snare operated by a multi-functional handle
- Minimizes procedure time by eliminating device exchanges between needle and snare
- Enables immediate snaring to avoid cushion/bleb absorption
- Allows for immediate injection for hemostasis - if needed after polyp resection
- Available in three unique snare options including the new multi-stage Lariat® snare



NEW
iSnare® system -
Lariat® snare



iSnare® system -
oval snare



iSnare® system -
hexagonal snare

For more information on the iSnare® system — contact US Endoscopy today.

product number	description	needle gauge	needle length	snare size	sheath diameter	sheath length	units/box
BX00711084	Lariat® snare	25	5mm	oval shape: 30 x 55mm hex shape: 10 x 28mm diamond shape: 6 x 15mm	2.5mm	230cm	5
BX00711085	oval snare	23	5mm	25 x 40mm	3.0mm	230cm	5
BX00711086	oval snare	25	5mm	25 x 40mm	3.0mm	230cm	5
BX00711088	hexagonal snare	25	5mm	25 x 40mm	3.0mm	230cm	5
BX00711089	hexagonal snare	23	5mm	25 x 40mm	3.0mm	230cm	5

Visit www.usendoscopy.com to see product videos, literature, and much more.

Distributed by:



listening... and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639-4494

fax +1 440 / 639-4495

email global@usendoscopy.com

www.usendoscopy.com