

iSnare[®] system



*Two snare styles
available*

Maintain scope positioning...

The iSnare[®] system from US Endoscopy is a family of unique, multi-functional devices ideal for endoscopic mucosal resection and saline assisted polypectomy.

The iSnare[®] system is designed to serve all of the following endoscopic needs:

- endoscopic mucosal resection
- saline assisted polypectomy
- hemostatic injection
- applying dye or stain during chromoendoscopy
- piecemeal resection of larger polyps
- tattooing pre-or post-polypectomy
- adherent clot removal

and any other procedures requiring injection therapy and/or snare resection techniques.

The iSnare[®] system...

- *Features alternating injection needle and polypectomy snare operated by a multi-functional handle*
- *Facilitates immediate injections in response to bleb absorption/dissipation or bleeding*
- *Minimizes wasted procedural time – inject or snare at any time throughout the procedure without device exchanges while maintaining scope position*
- *Features a spring-loaded, luer lock handle for consistent needle projection and retraction*
- *Offers protection against intraprocedural needle perforation of the catheter*
- *Allows the ability to work with many commonly used injection media*

Stiffer hexagonal snare style

For more information on the iSnare[®] system, contact US Endoscopy today.

product number	description	needle gauge	needle length	snare size	sheath diameter	sheath length	active cord connector	units/box
00711085	oval snare	23	5mm	2.5 x 4.0cm	3.0mm	230cm	Olympus	5
00711086	oval snare	25	5mm	2.5 x 4.0cm	3.0mm	230cm	Olympus	5
00711088	hexagonal snare	25	5mm	2.5 x 4.0cm	3.0mm	230cm	Olympus	5
00711089	hexagonal snare	23	5mm	2.5 x 4.0cm	3.0mm	230cm	Olympus	5

Visit www.usendoscopy.com to see product videos, literature, and much more.

Distributed by:



listening... and delivering solutions[®]

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639-4494

fax +1 440 / 639-4495

email global@usendoscopy.com

www.usendoscopy.com