

BioGuard®
air/water & suction valves



Improving infection prevention practices...

The single-use BioGuard® air/water & suction valves offer an alternative option to manual cleaning and reprocessing.



BX00711780
2-piece kit



BX00711783
3-piece kit



BX00711784
4-piece kit

The BioGuard® air/water & suction valves...

- Minimize the risk of cross-contamination between patients
- Provide consistent performance – use a new kit for each procedure
- Help eliminate tracking concerns in an infection outbreak*
- Support adherence to industry guidelines

For more information on the BioGuard® air/water & suction valves, contact US Endoscopy today.

product number	description	sterile	units/box
BX00711780	BioGuard® air/water & suction valves - 2 piece kit	yes	50
BX00711783	BioGuard® air/water & suction valves - 3 piece kit *kit includes BioShield® biopsy valve (BX00711129)	yes	50
BX00711784	BioGuard® air/water & suction valves - 4 piece kit *kit includes BioShield® biopsy valve (BX00711129) and Torrent® scope connector (BX00711570)	yes	50

Each kit individually packaged and sealed in tray.
Compatible with Olympus GI endoscopes (excluding EUS echoendoscopes that utilize a balloon channel).

Distributed by:



DIAGMED
HEALTHCARE

Making a difference to Health

Diagmed Healthcare Ltd, Lumley Close, Thirsk Industrial Park,
Thirsk, North Yorkshire, England, YO7 3TD
01845 526660
info@diagmed.co.uk
diagmed.healthcare



listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

*British Society of Gastroenterology Endoscopy Committee: BSG Guidance on For Decontamination of Equipment For Gastrointestinal Endoscopy. 3, 10, 50.
Olympus® is a registered trademark of Olympus Corporation.
US Endoscopy is a subsidiary of STERIS Corporation.
Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.
©2017 US Endoscopy. All Rights Reserved.