

bite blocks




A broad offering to fit all patient needs...

Our extensive line of bite blocks offers a variety of options in sizes, designs, and materials.


Standard

The sturdy, yet forgiving material of the standard bite blocks protects delicate instruments and provides comfort for patients.



maxi bite block / 60FR

- available with a velcro, latex-free or woven strap, or without a strap
- with dental retention rim



standard bite block / 48FR

- available with velcro and latex-free straps or without a strap
- with dental retention rim



intermediate bite block / 48FR

- latex-free strap
- without dental retention rim

product number	description	opening size	dental retention rim	units/box
00712801	standard - latex-free strap	16mm / 48 FR	yes	100
00712802	standard - no strap	16mm / 48 FR	yes	100
00712803	standard - velcro strap	16mm / 48 FR	yes	100
00712804	maxi - latex-free strap	20mm / 60 FR	yes	100
00712805	maxi - no strap	20mm / 60 FR	yes	100
00712806	maxi - velcro strap	20mm / 60 FR	yes	100
00712809	intermediate - latex-free strap	16mm / 48 FR	no	100
New 00712812	maxi - woven strap	20mm / 60 FR	yes	100

models to ensure the best for the patient and type of procedure.

Specialty

From a pediatric bite block to a purpose-built O₂ delivery bite block, our specialty bite blocks are designed to fit the unique needs of patients.



pediatric bite block / 38FR

- latex-free strap
- without dental retention rim



Vasto™ bite block / 72FR

- featuring larger side ports for scope repositioning and suction access
- latex-free strap
- with dental retention rim



Respa™ O₂ delivery bite block / 60 FR **New**

- can be used to house the nasal cannula during procedures to easily deliver O₂ to the patient orally
- latex-free strap
- with dental retention rim

For more information, contact US Endoscopy today.

product number	description	opening size	dental retention rim	units/box
00712808	pediatric - latex-free strap	12.6mm / 38 FR	no	25
00712810	Vasto™ bite block - latex-free strap	24mm / 72 FR	yes	100
New 00712811	Respa™ O ₂ delivery bite block	20mm / 60 FR	yes	100

Velcro® is a registered trademark of Velcro Industries B.V.

US Endoscopy is a subsidiary of STERIS Corporation. Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy



listening... and delivering solutions®

5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226
email global@usendoscopy.com

www.usendoscopy.com