


**US Endoscopy
injection needles**



A variety of needles for a variety of needs...

US Endoscopy features a full line of injection needles to suit nearly all of your injection needs - from saline assisted polypectomy and enteroscopy to esophageal varices, ulcers, and tattooing.

Carr-Locke injection needle – designed to operate in the most tortuous conditions

now available in 23 gauge

- Made with a stainless steel, spring sheath to virtually eliminate sheath kinking, while offering exceptional functionality and needle performance
- Ultra smooth distal metal hub is gentle on friable tissue, aids in tamponade, and prevents sheath perforations
- Spring-loaded handle automatically retracts the needle, preventing inadvertent needle sticks

product number	sheath diameter (mm)	length (cm)	needle projection (mm)	needle gauge	units/box
00711811	2.5	230	5	25	5
00711812	2.5	230	5	23	5
00711813 (Pentax compatible)	2.5	230	5	25	5

New

Carr-Locke injection needle

Articulator™ injection needle – ideal for everyday use and dependable for your most immediate needs

- Luer-lock, spring-loaded handle helps ensure the needle is either engaged or retracted, reducing the risk of inadvertent scope damage and incidental needle sticks
- Kink resistant, stainless steel sheath to pass smoothly and perform in fully articulated endoscopes

product number	sheath diameter (mm)	length (cm)	needle projection (mm)	needle gauge	units/box
00711803	2.5	160	4	25	5
00711804	2.5	160	5	25	5
00711807	2.5	230	5	25	5
00711808*	2.5	350	5	25	5

Articulator™ injection needle

Vari-Safe™ injection needle – cost effective and reliable, coupled with superior performance

- Luer-lock, spring-loaded handle helps promote maximum needle projection, while protecting against scope damage and inadvertent needle sticks

product number	sheath diameter (mm)	length (cm)	needle projection (mm)	needle gauge	units/box
00711819	2.5	230	5	23	1

Vari-Safe™ injection needle

Distributed by:



listening... and delivering solutions®
 5976 Heisley Road
 Mentor, OH 44060 USA
 phone +1 440 / 639-4494
 fax +1 440 / 639-4495
 email global@usendoscopy.com
www.usendoscopy.com

* 00711808 does not have the spring-loaded handle
 Listening... and delivering solutions®, the US Endoscopy listening... and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy
 PENTAX® is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha