

SmartBand® SafeGrip™ Multi-Band Ligation Kit

Scientifically designed to improve ligation performance




SmartBand® SafeGrip™ Multi-Band Ligation Kit

Features and Benefits

Anti-Slip Band Designed for Retention

SmartBand SafeGrip's Anti-Slip band design provides on average 33% more gripping capability,¹ to keep the band on the varix. The unique Anti-Slip band structure potentially leads to less band slippage to confidently treat the patient.

33% More Gripping Capability on Average¹



Predictable Band Deployment

The SmartBand handle incorporates human factors engineering for predictable band deployment.¹ With a natural twist of the wrist, each band is released with confidence. Tactile feedback upon band deployment provides an additional layer of assurance.



Features and Benefits (Continued)

Intuitive Setup

SmartBand improves on device setup with an intuitive hook and loop theme that helps to load the device correctly the first time.



The tear notches on the UV light protected pouch ensure easy access to the product.



SmartBand Multi-Band Ligation Pack

SmartBand Ligation Packs are available when more bands are needed to complete the procedure. The SmartBand Ligation Pack has the potential to reduce procedure time and cost by eliminating the need to open a new device.



Intended Use

Used to endoscopically ligate esophageal varices at or above the gastroesophageal junction and to ligate internal hemorrhoids.

Please refer to Instructions for Use included with each device.

Product Specifications

Order Number	Product Name	Scope Outer Diameter	Number of Bands	Band Material	Shelf Life (months)
SLK-6-LF	SmartBand SafeGrip Ligation Kit	8.6 – 11.6	6	Latex Free	24
SLP-6-LF	SmartBand SafeGrip Ligation Pack	8.6 – 11.6	6	Latex Free	24

¹ Data on file.

Diagmed Healthcare Ltd

Lumley Close, Thirsk Industrial Park, Thirsk,
North Yorkshire, UK, YO7 3TD

Tel: (01845) 526660 Fax: (01845) 522199
Email: info@diagmed.co.uk
Website: diagmed.healthcare



DIAGMED
HEALTHCARE