


Infinity[®] ERCP sampling device



*"Clinical study shows that the Infinity[®]
ERCP sampling device yields a
diagnostic sensitivity of up to 85%*."*


*T. Barrioz, M. Wangermez, P. Levillain, M. Beauchant. The Infinity[®] brush improves the results of brush cytology for malignant biliary strictures. Societe Nationale Francaise de Gastro-Enterologie (SNFGE). <http://www.snfge.org/jfhod2013/7407.html>[25-4-2013 10:13:01].



Looking to collect quality samples during ERCP?

The Infinity® ERCP sampling device is purpose built for collecting substantial and quality samples from strictures in the biliary duct.

The Infinity® device features...



Other features include...

- atraumatic tip
- nitinol drive wire
- salvage cytology
- long and short guidewire compatibility

For more information on the Infinity® ERCP sampling device, contact US Endoscopy today.

product number	sheath diameter (mm)	length (cm)	sterile	units/box
BX0071652	3 / 9FR	200	yes	5
BX0071653	2.5 / 7.5FR	200	yes	5

Distributed by:



listening...and delivering solutions®

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

760746 Rev. E

US Endoscopy is a subsidiary of STERIS Corporation.

© US Endoscopy or its affiliates. All rights reserved.

Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy

Follow us:    