



**ENDOSCOPIC TATTOO**  
*Take The Next Step In The Fight  
Against Colon Cancer*




**DIAGMED**  
HEALTHCARE




## WHAT IS SPOT® EX?

Spot® Ex is indicated for endoscopically tattooing lesions in the GI Tract for up to 36 months to facilitate both clinical surveillance and surgical localisation.



## PREPARATION



- 1 SHAKE IT
- 2 ATTACH IT
- 3 PRIME IT



“ Colonoscopic tattooing is performed to enable future identification, at colonoscopy or surgery, of malignant lesions (proven or suspected), polypectomy, EMR, or ESD sites, difficult-to-detect polyps, or dysplastic areas. **All such lesions**, other than those definitely located in the cecum, adjacent to the ileocecal valve, or in the low rectum, **should be tattooed.**”

A sterile and biocompatible pre-packaged suspension ... (Spot) has been developed...and this has enhanced the accessibility, ease of use, and safety of the procedure.”



Spot® Ex is much darker than Spot®, making it easier to find at the time of surgery

## SPOT® EX TATTOOS ARE SOCIETY RECOMMENDED:

## ORDER DETAILS



- **Product Code:** GIS-46 SPOT
- **Description:** Spot® Ex Endoscopic Tattoo
- **Box Size:** 10

Diagmed Healthcare Ltd  
Lumley Close, Thirsk Industrial Park, Thirsk  
North Yorkshire, UK, YO7 3TD

Tel: (01845) 526660 Fax: (01845) 522199  
Email: [info@diagmed.co.uk](mailto:info@diagmed.co.uk)  
Website: [diagmed.healthcare](http://diagmed.healthcare)

