

ESOPHAGEAL

1 ECN - Fully Covered Asymmetrical

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
ECN18090-X070	60	90	24-18-24	70	6 / 18	No
ECN18110-X070	80	110	24-18-24	70	6 / 18	No

2 ECV - Fully Covered Symmetrical

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
ECV20080-X070	40	80	26-20-26	70	6 / 18	No
ECV20110-X070	70	110	26-20-26	70	6 / 18	No
ECV20140-X070	100	140	26-20-26	70	6 / 18	No
ECV20170-X070	130	170	26-20-26	70	6 / 18	No

2 EPW - Partially Covered

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
EPW-20-080-X070	40	80	26-20-26	70	6 / 18	No
EPW-20-100-X070	60	100	26-20-26	70	6 / 18	No
EPW-20-120-X070	80	120	26-20-26	70	6 / 18	No
EPW-20-150-X070	110	150	26-20-26	70	6 / 18	No

2 EBN - Double Covered for Benign Indications

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
EBN22080-Z070	40	80	28-22-28	70	8 / 24	No
EBN22120-Z070	80	120	28-22-28	70	8 / 24	No
EBN24080-Z070	40	80	30-24-30	70	8 / 24	No
EBN24100-Z070	80	120	30-24-30	70	8 / 24	No

2 EPBA - Partially Covered TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
EPBA20070-C230	30	70	26-20-26	230	3.47 / 10.5	Yes
EPBA20090-C230	50	90	26-20-26	230	3.47 / 10.5	Yes
EPBA20110-C230	70	110	26-20-26	230	3.47 / 10.5	Yes
EPBA20130-C230	90	130	26-20-26	230	3.47 / 10.5	Yes
EPBA20150-C230	110	150	26-20-26	230	3.47 / 10.5	Yes

3 ECW - Fully Covered with Anti-Reflux Valve

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
ECW22090-Z070	50	90	28-22-28	70	8 / 24	No
ECW22120-Z070	80	120	28-22-28	70	8 / 24	No
ECW22160-Z070	120	160	28-22-28	70	8 / 24	No

3 EPV - Partially Covered with Anti-Reflux Valve

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
EPV22080-Z070	40	80	28-22-28	70	8 / 24	No
EPV22100-Z070	60	100	28-22-28	70	8 / 24	No
EPV22120-Z070	80	120	28-22-28	70	8 / 24	No
EPV22140-Z070	100	140	28-22-28	70	8 / 24	No

3 ECBB - Fully Covered for Bariatric Surgery Complications

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
ECBB-28-180-090	145	180	32-28-28	90	8 / 24	No
ECBB-28-210-090	175	210	32-28-28	90	8 / 24	No
ECBB-28-240-090	205	240	32-28-28	90	8 / 24	No

BILIARY

4 SHS - Uncovered FAST Delivery System **NEW**

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
SHS-08-040-S180	26	40	11.5-8-11.5	180	2.83 / 8.5	Yes
SHS-08-060-S180	46	60	11.5-8-11.5	180	2.83 / 8.5	Yes
SHS-08-080-S180	66	80	11.5-8-11.5	180	2.83 / 8.5	Yes
SHS-08-100-S180	86	100	11.5-8-11.5	180	2.83 / 8.5	Yes
SHS-10-040-S180	26	40	13.5-10-13.5	180	2.83 / 8.5	Yes
SHS-10-060-S180	46	60	13.5-10-13.5	180	2.83 / 8.5	Yes
SHS-10-080-S180	66	80	13.5-10-13.5	180	2.83 / 8.5	Yes
SHS-10-100-S180	86	100	13.5-10-13.5	180	2.83 / 8.5	Yes

4 SHS - Uncovered Endoscopic TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
SHS-08-040-180	26	40	11.5-8-11.5	180	2.33 / 7	Yes
SHS-08-060-180	46	60	11.5-8-11.5	180	2.33 / 7	Yes
SHS-08-080-180	66	80	11.5-8-11.5	180	2.33 / 7	Yes
SHS-08-100-180	86	100	11.5-8-11.5	180	2.33 / 7	Yes
SHS-10-040-180	26	40	13.5-10-13.5	180	2.33 / 7	Yes
SHS-10-060-180	46	60	13.5-10-13.5	180	2.33 / 7	Yes
SHS-10-080-180	66	80	13.5-10-13.5	180	2.33 / 7	Yes
SHS-10-100-180	86	100	13.5-10-13.5	180	2.33 / 7	Yes
SHS-10-120-180	106	120	13.5-10-13.5	180	2.33 / 7	Yes

4 SHS - Uncovered Percutaneous Delivery


	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
SHS-10-040-060	26	40	13.5-10-13.5	60	2.33 / 7	Yes
SHS-10-060-060	46	60	13.5-10-13.5	60	2.33 / 7	Yes
SHS-10-080-060	66	80	13.5-10-13.5	60	2.33 / 7	Yes
SHS-10-100-060	86	100	13.5-10-13.5	60	2.33 / 7	Yes
SHS-10-120-060	106	120	13.5-10-13.5	60	2.33 / 7	Yes

4 BNY - Uncovered Hilar Stenoses Endoscopic TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BNY-08-060-180	46	60	11.5-8-11.5	180	2.33 / 7	Yes
BNY-08-080-180	66	80	11.5-8-11.5	180	2.33 / 7	Yes
BNY-08-100-180	86	100	11.5-8-11.5	180	2.33 / 7	Yes
BNY-08-120-180	106	120	11.5-8-11.5	180	2.33 / 7	Yes
BNY-10-060-180	46	60	13.5-10-13.5	180	2.33 / 7	Yes
BNY-10-080-180	66	80	13.5-10-13.5	180	2.33 / 7	Yes
BNY-10-100-180	86	100	13.5-10-13.5	180	2.33 / 7	Yes
BNY-10-120-180	106	120	13.5-10-13.5	180	2.33 / 7	Yes

4 BNY - Uncovered Hilar Stenoses Percutaneous

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BNY-08-060-060	46	60	11.5-8-11.5	60	2.33 / 7	Yes
BNY-08-080-060	66	80	11.5-8-11.5	60	2.33 / 7	Yes
BNY-08-100-060	86	100	11.5-8-11.5	60	2.33 / 7	Yes
BNY-08-120-060	106	120	11.5-8-11.5	60	2.33 / 7	Yes
BNY-10-060-060	46	60	13.5-10-13.5	60	2.33 / 7	Yes
BNY-10-080-060	66	80	13.5-10-13.5	60	2.33 / 7	Yes
BNY-10-100-060	86	100	13.5-10-13.5	60	2.33 / 7	Yes
BNY-10-120-060	106	120	13.5-10-13.5	60	2.33 / 7	Yes



4 BCG - Fully Covered with Flaps FAST Delivery System **NEW**

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BCG-08-040-S180	20	40	8-8-8	180	2.83 / 8.5	Yes
BCG-08-060-S180	40	60	8-8-8	180	2.83 / 8.5	Yes
BCG-08-080-S180	60	80	8-8-8	180	2.83 / 8.5	Yes
BCG-10-040-S180	20	40	10-10-10	180	2.83 / 8.5	Yes
BCG-10-060-S180	40	60	10-10-10	180	2.83 / 8.5	Yes
BCG-10-080-S180	60	80	10-10-10	180	2.83 / 8.5	Yes

4 BCG - Fully Covered with Flaps Endoscopic TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BCG-10-040-180	20	40	10-10-10	180	2.83 / 8.5	Yes
BCG-10-060-180	40	60	10-10-10	180	2.83 / 8.5	Yes
BCG-10-080-180	60	80	10-10-10	180	2.83 / 8.5	Yes
BCG-10-100-180	80	100	10-10-10	180	2.83 / 8.5	Yes

4 BCG - Fully Covered with Flaps Percutaneous TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BCG-10-040-060	20	40	10-10-10	60	2.83 / 8.5	Yes
BCG-10-060-060	40	60	10-10-10	60	2.83 / 8.5	Yes
BCG-10-080-060	60	80	10-10-10	60	2.83 / 8.5	Yes
BCG-10-100-060	80	100	10-10-10	60	2.83 / 8.5	Yes

4 BCS - Fully Covered with Flaps and Flare Endoscopic TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BCS10040-D180	23	40	10-10-13.5	180	2.83 / 8.5	Yes
BCS10060-D180	43	60	10-10-13.5	180	2.83 / 8.5	Yes
BCS10080-D180	63	80	10-10-13.5	180	2.83 / 8.5	Yes
BCS10100-D180	83	100	10-10-13.5	180	2.83 / 8.5	Yes
BCS10120-D180	103	120	10-10-13.5	180	2.83 / 8.5	Yes

4 BCS - Fully Covered with Flaps and Flare Percutaneous TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BCS10040-R060	23	40	10-10-13.5	60	2.83 / 8.5	Yes
BCS10060-R060	43	60	10-10-13.5	60	2.83 / 8.5	Yes
BCS10080-R060	63	80	10-10-13.5	60	2.83 / 8.5	Yes
BCS10100-R060	83	100	10-10-13.5	60	2.83 / 8.5	Yes
BCS10120-R060	103	120	10-10-13.5	60	2.83 / 8.5	Yes

4 BCT - Fully Covered with Flaps and Flare and Lasso Endoscopic TTS Delivery

	Stent (mm)			Delivery Device		
	Usable length	Total length	Diameters	Length (cm)	Diameter (mm/ft)	Resheathable
BCT-10-040-180	23	40	10-10-13.5	180	3.4 / 10.2	Yes
BCT-10-060-180	43	60	10-10-13.5	180	3.4 / 10.2	Yes
BCT-10-080-180	63	80	10-10-13.5	180	3.4 / 10.2	Yes
BCT-10-100-180	83	100	10-10-13.5</			