

Falcon™ rotatable retrieval basket



Capture long objects from any angle...

The Falcon™ rotatable retrieval basket offers flexibility when retrieving long or oddly shaped objects such as toothbrushes, utensils, and pens.

The Falcon™ basket is...

- rotatable – enables the physician to capture objects easier
- single-use – eliminates reprocessing and ensures consistent performance
- three dimensional – allows clinician to retrieve objects at any angle

*“ASGE Guidelines for the management of ingested foreign bodies states, “Equipment that should be readily available includes forceps, polypectomy snare, polyp grasper, and **Dormier basket...**”*

(Gastrointestinal Endoscopy 6/2011)

4 wire design

For more information on the Falcon™ rotatable retrieval basket, contact US Endoscopy today.

product number	sheath diameter (mm)	length (cm)	basket size	units/box
00711141	2.5	230	3.5 x 6.0	5

Distributed by:

us endoscopy

*listening...and delivering solutions**

5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

760840 Rev. A