








essential foreign body equipment

	Roth Net® retriever - select	Roth Net® maxi retrieval net	Roth Net® mini retrieval net	Talon® grasping device	Raptor® grasping device	Guardus® esophageal/gastric overtube	Falcon® rotatable retrieval basket
							
	00711169/00711170 Net size - 3.0cm x 6.0cm	00711059/00711165 Net size - 4.0cm x 8.0cm	00711057/00711168 Net size - 2.0cm x 4.5cm	00711175 Length - 160cm	00711177/00711178 Length - 230cm/200cm	00711146/00711147 00711148/00711149 Length - 25cm/50cm	00711141 Basket size 3.5cm x 6.0cm
Larger blunt object removal (AA battery...)		✓				✓	
Smaller blunt object removal (buttons, coins...)	✓	✓	✓		✓	✓	
Long object removal (spoon, toothbrush...)						✓	✓
Sharp-pointed object removal (toothpicks, bones, paperclips...)					✓	✓	
Disk battery removal	✓		✓			✓	
Food bolus removal	✓			✓		✓	

are you prepared?

The ASGE Guideline for the management of ingested foreign bodies states, "Equipment that should be readily available includes **rat tooth and alligator forceps, polypectomy snare, polyp grasper, Dormier basket, retrieval net, overtubes of esophageal and gastric lengths, and a foreign body protector hood.**"

Gastrointestinal Endoscopy 6/2011

The above chart suggests applications for US Endoscopy foreign body management devices. Only persons having adequate training and proficiency in endoscopic foreign body removal technique are qualified to determine appropriate device usage in any given procedure.

product	qty.	complete foreign body inventory checklist														
Roth Net® retrieval nets																
Talon® grasping device																
Guardus® esophageal overtubes																
Guardus® gastric overtubes																
Falcon® retrieval basket																
Raptor® grasping devices																
polyp grasper																
foreign body protector hood																
initials & date																

product	qty.	complete foreign body inventory checklist														
Roth Net® retrieval nets																
Talon® grasping device																
Guardus® esophageal overtubes																
Guardus® gastric overtubes																
Falcon® retrieval basket																
Raptor® grasping devices																
polyp grasper																
foreign body protector hood																
initials & date																

US Endoscopy is a subsidiary of STERIS Corporation.
 Listening...and delivering solutions®, the US Endoscopy listening...and delivering solutions design, and all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.



phone (1) 440 / 639-4494
 fax (1) 440 / 639-4495
 customer service (1) 800 / 769-8226
www.usendoscopy.com