

# Endoscopy Irrigation Devices

To address cross-contamination concerns and support patient safety, Diagmed has a wide range of irrigation device solutions.



**DIAGMED**  
HEALTHCARE

MAKING A DIFFERENCE TO HEALTH



AQUASHIELD®  
CONNECTED TO  
OLYMPUS ELECTRIC PUMP



24-HOUR RED STICKER  
Which identifies the date and  
time of the setup

DUAL TUBE  
SYSTEM  
With easy pinch  
clamp

## AquaShield® CO2 Water Bottle System

The disposable AquaShield® system offers a safe and convenient alternative to reprocessing reusable water bottle systems.

### Features:

- Used in conjunction with the endoscope's air / water function to clean the scope lens
- Saves time and money associated with reprocessing – there is no need to reprocess or replace lost/damaged components
- Proactively supports patient safety by reducing the risk of bacterial contamination that may remain in reusable systems
- A 24-hour use system – use one for the day and throw it away
- Available for use with air or CO2 insufflation
- Dual tube system with easy pinch clamp
- Long shelf life Vs rival products

### Ordering Information

Product Code	Description	Compatibility	Box Size
00711541	AquaShield® CO2 Water Bottle System	Olympus 140 / 160 / 180 Series (10cm line)	25
00711543	AquaShield® CO2 Water Bottle System	Braun and Fresenius Compatible	25
00711550	AquaShield® Water Bottle System	Pentax	25
00711551	AquaShield® CO2 Water Bottle System	Pentax (152cm line)	25
00711549	AquaShield® CO2 Water Bottle System	Olympus (152cm line) & filter	25
00711553	AquaShield® CO2 Water Bottle System	Fujifilm	25
00711599	AquaShield® Holder	-	1



## Disposable Washing Catheter

Our Disposable Washing Catheter is constructed of high quality PTFE Teflon with a luer-lock connector. PTFE Teflon has one of the lowest coefficients of friction of any copolymer on the market today. This enables the Endoscopist to pass this catheter down the biopsy channel of an endoscope with very little resistance.

### Features:

- Offers superior flushing ability compared to some other methods
- Suitable for use with all endoscopes with an accessory channel diameter of at least 2.8mm
- Easy to use and fully disposable

### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Catheter Length (cm)	Box Size
WC-1	Disposable Washing Catheter	2.3	240	1

## Torrent® Irrigation Tubing

The Torrent® Irrigation Tubing is compatible with Olympus, ERBE, and Byrne electric pumps.

### Features:

- Offers a single-use connector, supporting patient safety and reducing risk of cross-contamination
- Features 24-hour use irrigation tubing – use for a day then throw it away
- Eliminates time and costs associated with reprocessing and replacing reusable tubing

### Ordering Information


Product Code	Description	Compatibility	Length (cm)	Box Size
00711562	Torrent® Irrigaton Tubing	For Olympus Electric Pumps	366	25
00711563	Torrent® Irrigaton Tubing	For ERBE / Byrne Electric Pumps	366	25



**SINGLE USE CONNECTOR**  
Reducing the risk of cross-contamination



**TORRENT® SCOPE CONNECTOR FOR PENTAX SCOPES**



## EndoShield CO2 Water Bottle System

Fully disposable 24-hour water bottle caps with CO2 compatibility. Quick, easy, reliable and great value-for-money!



CONTROL VALVE ON CO2 CONNECTION PORT



### Features:

- Control valve on CO2 connection port which prevents back flow of water to CO2 insufflator. It also prevents air from escaping when using room air instead of CO2
- Features an easy pinch clamp which reduces water leaks and drips
- Compatible with a large range of water bottles
- 24-hour date and time sticker
- Greater scope compatibility which reduces the need to reprocess and replace expensive reusable systems
- Longer 36 month shelf life
- Fully disposable product which reduces the risk of cross contamination and delivers sterile water to the scope safely and efficiently

### Ordering Information

Product Code	Description	Compatibility	Box Size
ENDOSHEILDFB	EndoShield CO2 Water Bottle System	For Fujinon Scopes and Baxter Bottles	10
ENDOSHIELDFF	EndoShield CO2 Water Bottle System	For Fujinon Scopes and Fresenius Bottles	10
ENDOSHIELDOB	EndoShield CO2 Water Bottle System	For Olympus Scopes and Baxter Bottles	10
ENDOSHIELDFOF	EndoShield CO2 Water Bottle System	For Olympus Scopes and Fresenius Bottles	10

## EndoFlow 24 Hour Irrigation Tubing

Fully disposable 24-hour irrigation tubing. Quick, easy, reliable and great value-for-money!

### Features:

- Fully disposable product which reduces the risk of cross-contamination and delivers sterile water through the scope safely and efficiently
- Features a 24-hour date / time label
- Greater scope compatibility for Olympus and Fujinon
- Longer (3m) tubing which caters for most procedure room configurations

### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Box Size
ENDOFLOW-B	24-Hour Irrigation Tubing	Baxter Bottles	10
ENDOFLOW-F	24-Hour Irrigation Tubing	Fresenius Bottles	10
ENDOFLOW ADAPTOR	Biopsy Port Channel Adaptor	Universal Fit	10

## EndoFlow CO2 Insufflator Tubing Set

The new EndoFlow CO2 reduces the risk of cross-contamination with our disposable CO2 tubing system.

### Features:

- Fully disposable product which reduces the risk of cross-contamination and delivery CO2 via the scope safely and efficiently
- Compatible with many leading CO2 insufflators
- Compatibility with any disposable water bottle system making it extremely versatile
- 24-hour date and time sticker
- Longer 36 month shelf life

### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Box Size
ENDOFLOW CO2	EndoFlow CO2 Insufflator Tubing Set	Universal	1



ENDOFLOW CO2  
CO2 Insufflator Tubing Set



ENDOFLOW-F  
24-Hour Irrigation Tubing



Lumley Close, Thirsk Industrial Park, Thirsk, North Yorkshire, England, YO7 3TD  
01845 526660 | [info@diagmed.co.uk](mailto:info@diagmed.co.uk) | [diagmed.healthcare](http://diagmed.healthcare)

