


Lariat®
snare




Address various polyp sizes with one snare...


The multifunctional design of the Lariat® snare offers hot/cold capability to address a variety of resection needs, including piecemeal or en bloc resection of multiple polyps.

The Lariat® snare features...


- Three distinct sizes and shapes eliminating the need to exchange snares
- Snare loop design maintains its shape and integrity after multiple resections
- Stiffer snare wire minimizes distal deflection



Stage 1:
30mm oval



Stage 2:
10mm hex



Stage 3:
6mm diamond

For more information on the Lariat® snare, contact US Endoscopy today.

product number	description	sheath diameter (mm)	length (cm)	approximate snare size (width x length)	units/box
BX00711119	Lariat® snare	2.4	230	oval shape: 30x55mm hex shape: 10x28mm diamond shape: 6x15mm	10
BX00711096	Lariat® snare	2.4	230	oval shape: 30x55mm hex shape: 10x28mm diamond shape: 6x15mm	20
BX00711097	Lariat® snare	2.4	230	oval shape: 30x55mm hex shape: 10x28mm diamond shape: 6x15mm	50

Distributed by:



5976 Heisley Road
Mentor, OH 44060 USA

phone +1 440 / 639.4494

fax +1 440 / 639.4495

email global@usendoscopy.com

www.usendoscopy.com

US Endoscopy is a subsidiary of STERIS Corporation.
All marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or trademarks owned by US Endoscopy.
© 2018 US Endoscopy or its affiliates. All rights reserved.

Follow us:    