

Endoscopic
Set

made in
Germany

HaemoCer™ Endoscopic SET

Haemostatic powder for
endoscopic application

- ✓ Tube length 2400 mm
- ✓ Outer tube diameter of 2.5 mm
- ✓ HEPA filter -
High Efficiency Particulate Air
- ✓ Usable instrument channel:
2.8 mm and above



DIAGMED
HEALTHCARE

BioCer

Endoscopic SET

Haemostatic Powder + Endoscopic Applicator



Haemostatic Powder

- Absorbable Polysaccharide Haemostat
- Containing no animal nor human components
- Plant-based polymer crosslinking that creates ultra-hydrophilic, biocompatible particles.
- Rapidly resorbed by amylase within 48h without leaving any residues behind
- Rapid off-the-shelf deployment, no special storage, defrosting or premixing



HaemoCer™ Endoscopic Applicator

- For upper and lower gastrointestinal bleedings
- Even effective in difficult to access regions
- No blocking in case of contact with moisture
- Tube length 2400 mm
- Outer tube diameter of 2.5 mm
- HEPA filter - High Efficiency Particulate Air
- Usable instrument channel: 2.8 mm and above

Ordering Information

Product Code	Description	Box Size
HAEMOCER ENDOSET	Endoscopic Applicator + HaemoCer™ PLUS 3g	1

🏠 **DIAGMED HEALTHCARE LTD**
Lumley Close, Thirsk Industrial Park, Thirsk
North Yorkshire, UK, YO7 3TD

☎ Tel: (01845) 526660 📠 Fax: (01845) 522199

✉ Email: info@diagmed.co.uk

🌐 Website: diagmed.healthcare



DIAGMED
HEALTHCARE