

iSnare[®] system



“ The advantage of the combined snare over the braided snare, as seen in our study, is likely multifactorial. This study provides evidence that the integrated needle-snare may be superior to the snare alone for the removal of large, flat polyps.* ”

*T. Woodward, J. Crook, M. Raimondo, M. Wallace. "Improving complete EMR of colorectal neoplasia: a randomized trial comparing snares and injectate in the resection of large sessile colon polyps." GIE Journal online, 2015.

 **STERIS**


Endoscopy

Combining a needle and snare into one device to simplify advanced polyp resection...

The iSnare® system is a family of unique, multi-functional devices ideal for minimizing procedure time and optimizing clinical outcomes during submucosal injection polypectomy procedures including EMR (endoscopic mucosal resection) and SAP (saline assisted polypectomy).

The iSnare system

- » Features an alternating injection needle and snare operated by a multi-functional handle
- » Minimizes procedure time by reducing device exchanges between needle and snare
- » Enables immediate snaring to avoid cushion/bleb absorption
- » Allows for immediate injection for hemostasis - if needed after polyp resection
- » Available in three unique snare options including the multi-stage Lariat snare



For more information on the iSnare system — contact STERIS Endoscopy today.

product number	description	needle gauge	needle length (mm)	snare size (mm)	sheath diameter (mm)	sheath length (cm)	units/box
BX00711084	Lariat snare	25	5	oval shape: 30 hex shape: 10 diamond shape: 6	2.5	230	5
BX00711085	oval snare	23	5	25	3.0	230	5
BX00711086	oval snare	25	5	25	3.0	230	5
BX00711088	hexagonal snare	25	5	25	3.0	230	5
BX00711089	hexagonal snare	23	5	25	3.0	230	5

Visit www.steris.com to see product videos, literature, and much more.

Distributed by:



©2020 STERIS Endoscopy or its affiliates. All rights reserved.
Unless otherwise indicated, all marks denoted with © or ™ are registered with the U.S. Patent and Trademark Office, or are owned by STERIS or its affiliates.
Olympus® is a registered trademark of Olympus Optical, Inc.

STERIS
5960 Heisley Road
Mentor, OH 44060-1834 – USA
Customer Service: 800-548-4873
Fax: 440-639-4495
www.steris.com

Follow us: