



DIAGMED  
HEALTHCARE



# Roth Net®

## Speciality Retrieval

MAKING A DIFFERENCE TO HEALTH

# Roth Net® Retriever Portfolio

Roth Net® has been a trusted partner in foreign body and polyp retrieval for 20 years. Our broad portfolio of 20 different Roth Net® retrievers consists of sterile and non-sterile devices in a variety of sizes and shapes supporting polyp retrieval to the toughest foreign body and food bolus cases.

Device	NEW	Roth Net® 360 Retriever	Roth Net® Standard Retriever	Roth Net® Platinum™ Food Bolus Retriever	Roth Net® Platinum™ Universal Retriever	Roth Net® Select Retriever	Roth Net® Polyp Retriever	Roth Net® Mini Polyp Retriever	Roth Net® Maxi Retriever	Roth Net® Mini Retriever	Roth Net® Enteroscopic Retriever
	00711197	00711050/1052/1053/1162	00711155	00715050/5060	00711169/1170	00711150/1151	00711158/1159	00711059/1165	00711057/1168	00711152	
Objects											
SIZE	FOREIGN BODY RETRIEVAL					POLYP RETRIEVAL		SPECIALITY RETRIEVAL			
	Foreign Bodies / Food Bolus / Polyps	Round / Blunt Foreign Bodies	Food Bolus	Most Foreign Bodies	Food Bolus / Most Foreign Bodies	Polyps	Polyps	Large Foreign Bodies / EMR Specimens	Small Foreign Bodies	Foreign Bodies / Polyps in small bowel	
FEATURES	<ul style="list-style-type: none"> <li>360 degree controlled rotation</li> <li>Reinforced net offers strength and durability for tough foreign bodies</li> <li>Fully functions in retroflex or tortuous positions</li> </ul>	<ul style="list-style-type: none"> <li>Strong / reinforced net</li> </ul>	<ul style="list-style-type: none"> <li>Large / octagonal shape</li> <li>Opens to maximum capacity over multiple passes</li> </ul>	<ul style="list-style-type: none"> <li>Large / octagonal shape</li> <li>Strong / reinforced net for durability</li> </ul>	<ul style="list-style-type: none"> <li>Alternative net-weave for enhanced durability</li> <li>Unique netting for improved visibility</li> </ul>	<ul style="list-style-type: none"> <li>More transparent netting</li> <li>Easily encapsulates excised polyps</li> </ul>	<ul style="list-style-type: none"> <li>Round shape with memory points for consistent shape after multiple passes</li> </ul>	<ul style="list-style-type: none"> <li>Largest capacity retrieval net on the market</li> <li>Specifically designed for large foreign bodies and polyps</li> </ul>	<ul style="list-style-type: none"> <li>Designed for smaller foreign bodies</li> <li>Compatible with smaller channel endoscopes</li> </ul>	<ul style="list-style-type: none"> <li>For small bowel retrieval</li> </ul>	

# SETTING THE STANDARD

## Trusted Security

Our ROTH NET® retriever portfolio is proud to celebrate 25 years of evolution. As a pioneer, ROTH NET® is the very first retrieval net to market and has since led the innovation of foreign body and polyp retrievers.

ROTH NET® retrieval nets offer a durable solution for challenging retrieval procedures, including food bolus, foreign body and tissue retrieval. Today, our ever-evolving portfolio of 20 ROTH NET® retrievers continue to be the **Gold Standard** in endoscopic retrieval.

## ROTH NET® Select Polyp Retriever

The ROTH NET® Select Polyp Retriever is specifically designed to retrieve excised polyps, tissue samples and foreign bodies during flexible and rigid endoscopy procedures.

### Features:

- Maintains its shape and functionality throughout the procedure
- When wet, the ROTH NET® develops a slightly sticky feel to increase retention of resected polyp fragments when opening and closing during multiple sample collection
- Alternative weave design offers enhanced durability and strength for both foreign body and food bolus retrieval
- Unique netting for improved visibility throughout the procedure
- Available in sterile or non-sterile options

### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00711169	ROTH NET® Select Polyp Retriever - Sterile	2.5	230	3 x 6	5
00711170	ROTH NET® Select Polyp Retriever	2.5	230	3 x 6	5



*“Netting is tied to the wire in multiple locations, ensuring the net does not rupture or deform”*

*“ROTH NET® does not injure the mucosa and the net remains sturdy throughout the procedure”*



**UNIQUE NETTING**  
For improved visibility throughout the procedure



# ROTH NET® Retriever 360

The ROTH NET® Retriever 360 is the newest innovation and addition to the extensive ROTH NET® Retriever portfolio - offering best-in-class features. Featuring 360 degree, one-to-one rotation and a new spade shaped snare form that opens to full capacity even after multiple passes.




### Features:

- New spade shaped snare form that opens to full capacity even after multiple passes
- An alternative net weave offering durability, strength and visibility throughout the procedure
- 360 degree, one-to-one rotation allowing for controlled net placement and alignment


### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00711197	ROTH NET® Retriever - 360	2.5	230	3 x 5	5
00711198	ROTH NET® Retriever - 360 (Sterile)	2.5	230	3 x 5	10

360 DEGREE  
One-to-one rotation  
allowing for net  
placement and alignment



NEW SPADE SHAPE  
That opens to full capacity  
even after multiple passes



# ROTH NET® Standard Foreign Body Retriever

ROTH NET® Standard Foreign Body Retrievers are engineered with a strong, reinforced net material for extra security when passing the cricopharyngeus.

### Features:

- Several specialised designs offer optimal retrieval Nets for unique situations
- Patented Net design delivers a flexible pouch, ideal for retrieval within the GI Tract
- ROTH NET® retrievers are ideal for the removal of round, blunt or otherwise hard to grasp foreign bodies, as well as polyps and food bolus fragments




ROTH NET®  
Leading the way for 25 years!

*“A more durable retrieval net able to perform across multiple passes during challenging procedures and improved visibility when retrieving targeted objects”*

### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00711050	ROTH NET® Standard Foreign Body Retriever	2.5	230	3 x 6	5
00711052	ROTH NET® Standard Foreign Body Retriever	2.5	230	3 x 6	10
00711053	ROTH NET® Standard Foreign Body Retriever	2.5	160	3 x 6	5
00711162	ROTH NET® Standard Foreign Body Retriever Sterile	2.5	230	3 x 6	10
00711150	ROTH NET® Standard Polyp Retriever	2.5	230	3 x 6	5
00711151	ROTH NET® Standard Polyp Retriever	2.5	230	3 x 6	10



ROTH NET® MINI POLYP RETRIEVER

ROTH NET® MINI PAEDIATRIC RETRIEVER

## ROTH NET® Mini Polyp Retriever

The ROTH NET® Mini Polyp Retriever features a smaller net design, offering enhanced maneuverability and visualisation while retrieving polyps within the small confines of the colon.

### Features:

- Features a rounder snare shape that improves maneuverability and allows for a clear view during retrieval
- Contains memory points in the snare that enable the net to open to a full capacity and maintain a consistent shape during multiple uses

### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00711158	ROTH NET® Mini Polyp Retriever	2.8	230	3.5 x 2.5	10
00711159	ROTH NET® Mini Polyp Retriever Sterile	2.8	230	3.5 x 2.5	5

## ROTH NET® Mini Paediatric Foreign Body Retriever

The ROTH NET® Mini Paediatric Foreign Body Retriever is extremely effective for removing coins, watch batteries and a variety of other foreign bodies children tend to swallow which need urgent removal.

### Features:

- For use with gastroscopes that have a minimum 2.0mm working channel

### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00711057	ROTH NET® Mini Paediatric Foreign Body Retriever	1.8	160	2 x 4.5	5



“Versatile sample collection”

## ROTH NET® Platinum™ Food Bolus Retriever

The ROTH NET® Platinum™ Food Bolus Retriever is specifically designed for the challenges encountered in food bolus removal. Its octagonal net is designed to open to maximum capacity while maintaining its form within the Oesophagus.

### Features:

- Featuring larger capacity, ergonomic design and superior performance
- Unique octagonal design to maintain shape and open to maximum capacity
- Spring-like action of the ROTH NET® Platinum™ series offers a retrieval net with a sturdier form and shape across multiple passes



Upper edge of the food impaction  
Net partially opening while retracting the tip of the sheath



Encircling greater than 50% of the food bolus  
ROTH NET® and encased food bolus were slowly and completely withdrawn

“The food bolus was partially grasped and dislodged with the retrieval net. The net was then maneuvered to encircle the impaction. The food bolus was removed without any complication”

Professor Bennett E Roth M.D,  
David Geffen School of Medicine, UCLA


### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00711155	ROTH NET® Platinum™ Food Bolus Retriever	2.5	160	4 x 5.5	5





ROTH NET® MAXI FOREIGN BODY RETRIEVER



ROTH NET® ENTEROSCOPIC FOREIGN BODY RETRIEVER

## ROTH NET® Maxi Foreign Body Retriever

Ideal for large foreign body removal, such as batteries and EMR specimens too large for any other sized retriever or polyp trap.

### Features:

- For use in gastroscopes with 3.2 mm accessory channel or larger

*“Ability to utilise both sides of the ROTH NET® to aid with sample or foreign body collection”*

## ROTH NET® Enteroscopic Foreign Body Retriever

Ideal retrieval device for removal of excised polyps and foreign bodies, its unique net design offers a comfortable measure of safety and security during procedures. Extra-long 350cm catheter for use with enteroscopes.

### Ordering Information


Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00711059	ROTH NET® Maxi Foreign Body Retriever	3	160	4 x 8	5
00711165	ROTH NET® Maxi Foreign Body Retriever Sterile	3	230	4 x 8	5
00711152	ROTH NET® Enteroscopic Foreign Body Retriever	2.5	350	3 x 6	5

## ROTH NET® Platinum™ Polyp Retriever

The ROTH NET® Platinum™ Polyp Retriever features a unique spring-like design, as well as more transparent mesh for easy visibility during retrieval.

### Features:

- An oversized retrieval pouch which offers increased net capacity for large and piecemeal resected polyps, as well as multiple polyp collection
- For use in scopes with a minimum 3.2mm working channel
- The spring-like memory points of the ROTH NET® Platinum™ series offers a retrieval net with a sturdier form which maintains its shape during multiple openings in lengthy procedures



ROTH NET® PLATINUM™ POLYP RETRIEVER


## ROTH NET® Platinum™ Universal Retriever

Capable of addressing a variety of retrieval needs, this robust retriever is an all-encompassing solution for foreign bodies, food impactions and polyp retrieval.

### Features:

- Larger capacity, ergonomic design and superior performance
- The spring-like memory points of the ROTH NET® Platinum™ series offers a retrieval net with a sturdier form which maintains its shape during multiple openings in lengthy procedures

ROTH NET® PLATINUM™ UNIVERSAL RETRIEVER



### Ordering Information

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
00715150	ROTH NET® Platinum™ Polyp Retriever	3	230	4 x 5.5	5
00715160	ROTH NET® Platinum™ Polyp Retriever	3	230	4 x 5.5	5
00715050	ROTH NET® Platinum™ Polyp Retriever	2.5	230	4 x 5.5	5
00715060	ROTH NET® Platinum™ Universal Polyp Retriever	2.5	230	4 x 5.5	5

# Ordering Information for the ROTH NET® Portfolio

Product Code	Description	Catheter Diameter (mm)	Length (cm)	Net Size (cm)	Box Size
<i>Foreign Body Retrieval</i>					
00711197	ROTH NET® 360 Retriever - <b>NEW</b>	2.5	230	3 x 5	5
00711198	ROTH NET® 360 Retriever (Sterile) - <b>NEW</b>	2.5	230	3 x 5	10
00711050	ROTH NET® Standard Foreign Body Retriever	2.5	230	3 x 6	5
00711052	ROTH NET® Standard Foreign Body Retriever	2.5	230	3 x 6	10
00711053	ROTH NET® Standard Foreign Body Retriever	2.5	160	3 x 6	5
00711162	ROTH NET® Standard Foreign Body Retriever Sterile	2.5	230	3 x 6	10
00711155	ROTH NET® Platinum™ Food Bolus Retriever	2.5	160	4 x 5.5	5
00715050	ROTH NET® Platinum™ Polyp Retriever	2.5	230	4 x 5.5	5
00715060	ROTH NET® Platinum™ Universal Polyp Retriever	2.5	230	4 x 5.5	5
00711169	ROTH NET® Select Polyp Retriever Sterile	2.5	230	3 x 6	5
00711170	ROTH NET® Select Polyp Retriever	2.5	230	3 x 6	5
<i>Tissue Retrieval</i>					
00711150	ROTH NET® Standard Polyp Retriever	2.5	230	3 x 6	5
00711151	ROTH NET® Standard Polyp Retriever	2.5	230	3 x 6	5
00715150	ROTH NET® Platinum™ Polyp Retriever	3	230	4 x 5.5	5
00715160	ROTH NET® Platinum™ Polyp Retriever	3	230	4 x 5.5	5
<i>Speciality Retrieval</i>					
00711057	ROTH NET® Mini Paediatric Foreign Body Retriever	1.8	160	2 x 4.5	5
00711059	ROTH NET® Maxi Foreign Body Retriever	3	160	4 x 8	5
00711152	ROTH NET® Endoscopic Foreign Body Retriever	2.5	230	3 x 6	5
00711158	ROTH NET® Mini Polyp Retriever	2.8	230	3.5 x 2.5	10
00711159	ROTH NET® Mini Polyp Retriever Sterile	2.8	230	3.5 x 2.5	5
00711165	ROTH NET® Maxi Foreign Body Retriever Sterile	3	230	4 X 8	5



Lumley Close, Thirsk Industrial Park, Thirsk, North Yorkshire, England, YO7 3TD  
01845 526660 | info@diagmed.co.uk | diagmed.healthcare

