

Endoscopy Accessories

Protect patients and preserve
equipment with our range of useful
Endoscopic Ancillary Devices



DIAGMED
HEALTHCARE

a STERIS company

MAKING A DIFFERENCE TO HEALTH



**REINFORCED HAND
HOLDS:**
Allows easier transfer

Slyder™ Patient Transfer Sheet

Single use patient transfer device can slide up to 424Kgs (dry) / 415Kgs (wet).

Features and Benefits:

- Reinforced hand holds for better grip and reducing back strain
- Dimensions: 195cm X 85cm (standard)
- Latex free
- Reduces risk of injury if used properly
- Fully disposable for minimal contamination risk and requires no clean up



“Suitable for use in Theatres, GI Suite, Maternity, Radiology or Paediatrics”

Ordering Information

Product Code	Description	Dimensions	Box Size
SLYDER	Slyder™ Patient Transfer Sheet	195cm x 85cm	50 Units

Scope Tote - Infection Control Endoscope Transport Bag

The Scope Tote Bag is a transparent 24" X 20" re-sealable bag that easily holds all types of endoscopes as well as other types of soiled medical equipment - neatly, clearly and compactly.



SCOPE TOTE BAG

Features:

- Unlike a tray, Scope Tote bags completely contain the post-procedural equipment, reducing the risk of spreading infection
- Transparent for easy identification of contents while not exposing you to contaminants
- Each bag incorporates label for traceability and identification

Ordering Information

Product Code	Description	Size	Quantity
GIS-16	Scope Tote Bags	24" W X 20" H	25

Boost™ Head Positioner

The Boost™ Head Positioner is designed to help stabilise a patient's head during ERCP and other lengthy procedures, cradling it to enable easy access to the mouth for scope insertion and control throughout the procedure.



Features:

- Promotes and helps proper patient alignment and comfort during long procedures
- Cushions the patient's head which helps avoid pressure points
- Works for patient's in prone, semi-prone, and left-lateral positions
- Minimises cleaning out time and costs
- Minimises the risk of cross contamination and infection

BOOST™ HEAD POSITIONER



Ordering Information

Product Code	Description	Size	Box Size
00711654	Boost™ Head Positioner	H: 9.5cm X W: 19cm X L: 123cm	5

Endo-Boot™

Help protect your endoscopes from potentially costly repairs and downtime. Endo-boot™ protects the bending section and distal lenses of the endoscope during transportation and storage.

Features:

- Features a simple, easy to use design
- Reduce the chance of damage and costly repairs to the scope
- Disposable and suited to different sizes of endoscopes



ENDO-BOOT™ ON ENDOSCOPE

Ordering Information

Product Code	Description	Working Length	Pouch Dimensions	Opening Dimensions	Box Size
00711621	Endo-Boot™ Endoscope Tip Protector for Small Endoscopes	14cm	4cm	3mm-7mm	50
00711622	Endo-Boot™ Endoscope Tip Protector for Large Endoscopes	14cm	4cm	7mm-18mm	50


Keeper™ Device Clip

Stay organised and safe during procedures with the Keeper™ Device Clip.

Features:

- Keeps device neatly coiled
- Secures device catheter in between uses
- Protects and secures devices from kinking during procedures
- Can reduce the risk of cross-contamination by preventing accidental contact
- Ensures safe storage between users

Can be used with the Mio® organiser



KEEPER™ CLIP SECURING CATHETER

Ordering Information

Product Code	Description	Box Size
00711751	Keeper™ Device Clip	100

MIO® - Medical Device Organiser

Staying organised just got easier. The Mio® Medical device organiser is ideal for storing medical devices and supplies that can inadvertently fall on the floor, get misplaced, tangled or damaged during procedures.

Features:

- Easy setup - keeping devices organised and ready to go in one location
- Includes three oversized, transparent pouches to house a variety of accessories
- Wide pouches minimise wear and tear and prevent kinking from coiling devices too tightly
- Tape to fix towel on trolley, bedside rail, counter top or under mattress as needed
- Can reduce the risk of cross-contamination



“Device access and organisation made easy when you most need it”

MIO® ON COUNTER



Ordering Information

Product Code	Description	Dimensions (cm)	Pouch Dimensions (cm)	Box Size
00711750	Mio® Medical Device Organiser	86 X 76	25 X 28	20
TOWELS	ERCP Instrument Storage Towel	90 X 75	20 X 20	100



Lumley Close, Thirsk Industrial Park, Thirsk, North Yorkshire, England, YO7 3TD
01845 526660 | info@diagmed.co.uk | diagmed.healthcare

