

Endoscopy Garments

Extra comfort for your procedures, extra
care for your patient's peace of mind



DIAGMED
HEALTHCARE

a STERIS company

MAKING A DIFFERENCE TO HEALTH

D-Gowns - Paper Gown



Features:

- Soft and comfortable gown providing protection from cleaning fluid spray and spillages
- Open back provides ventilation
- The cool paper-type material is proven to reduce the problems associated with heat irritation to skin
- Water resistant
- Full arm-length cover with elasticated cuffs
- Two sizes: Large and Extra Large

Ordering Information

Product Code	Description	Box Size
D-GOWNS-BD	Paper Gowns	100
D-GOWNS-XL	Paper Gowns (Extra Large)	100

Ideal for screening programme

Dignity Colopants

Features:

- Dignity shorts specifically designed to maintain patient modesty before, during and after a lower GI procedure
- Opening controlled by Velcro fastening
- Individually packed
- Two sizes available - Large and Extra Large
- Unisex style
- Single-use and fully disposable to reduce the risk of cross contamination and infection
- Soft elasticated waist designed with patient comfort in mind

UNISEX
Dignity Colopants



Ordering Information

Product Code	Description	Size	Box Size
COLOPANTS-L-BD	Disposable Colonoscopy Pants	Large	100
COLOPANTS-XL-BD	Disposable Colonoscopy Pants	Extra Large	100

GIS-13 Long Moonpants

Features:

- Preserves patient's dignity and modesty – completely opaque soft material
- The most complete cover for patients with full leg cover
- Ideal for colonoscopy, sigmoidoscopy, barium enema and prostate ultrasound
- Adhesive tabs provide easy and safe opening and closing
- Single use and fully disposable



ADHESIVE TABS

Ordering Information

Product Code	Description	Available Sizes	Box Size
GIS-13	Long Moonpants	Medium / Large	50

PATIENT BIB - Disposable Endoscopy Bib

Features:

- Water resistant patient bib for any Upper GI endoscopy procedure
- Suitable for adult and paediatric procedures
- Includes large spill-resistant pocket for the collection of secretion, spittle and other bodily fluid
- Protects patient clothing and reduces spill risk for clinical staff
- Reduces the risk of cross contamination and infection



PATIENT BIB
Large pocket

Ordering Information

Product Code	Description	Box Size
PATIENT-BIB-BD	Disposable Endoscopy Bib	100

CONSTRUCTED WITH A SOFT, COMFORTABLE FABRIC

CHOICE OF LATEX-FREE ELASTIC WRIST OR SOFT CUFF



Defenz® Gown with Soft Cuff



Features:

- Soft, comfortable fabric with elasticated cuff
- Protects healthcare professionals from hazardous fluids
- The apron is designed to be easy to put on and easy to remove
- Protects the entire front of body and arms whilst the open back helps to keep the wearer cool and comfortable
- Water resistant
- Single use and fully disposable
- Premium PPE for the clean room

“Simply a softer approach.”

Ordering Information

Product Code	Description	Sizes Available	Box Size
GIS-8MLC-S	Defenz® Gown With Soft Cuff	Medium / Large	50

XGuard™ Radiation Reducing Gloves

Radiation reducing gloves have the same anatomic shape, size, elasticity and comfort typical of surgical gloves.

Features:

- Lead oxide impregnated to reduce the amount of ionising radiation reaching the hands by up to 55%
- Only 0.22mm or 0.30mm in thickness for optimal dexterity
- Latex deprived of allergic proteins minimise the risk of allergies
- Synthetic interior layer to further reduce allergic reactions
- Wide range of sizes from 6-9 including half sizes for a better fit
- Ideal where the user is regularly exposed to ionising radiations that may be encountered in the following types of procedures:
 - » Cardiology
 - » Nuclear Medicine
 - » Endoscopy
 - » Radiology
 - » Neuroradiology
 - » Orthopaedics
 - » Haemodynamic
 - » Resuscitation
 - » Pain Management
 - » Intensive Care
 - » Urology
 - » EP/Cath Lab

ELASTICITY AND COMFORT



Ordering Information

Product Code	Description	Size	Thickness	Box Size
RR - 1/6	Radiation Reducing Gloves	Size 6 (Small)	0.22mm	5
RR - 1/6.5	Radiation Reducing Gloves	Size 6.5 (Small)	0.22mm	5
RR - 1/7	Radiation Reducing Gloves	Size 7 (Medium)	0.22mm	5
RR - 1/7.5	Radiation Reducing Gloves	Size 7.5 (Medium)	0.22mm	5
RR - 1/8	Radiation Reducing Gloves	Size 8 (Large)	0.22mm	5
RR - 1/8.5	Radiation Reducing Gloves	Size 8.5 (Large)	0.22mm	5
RR - 1/9	Radiation Reducing Gloves	Size 9 (X-Large)	0.22mm	5
RR - 2/6	Radiation Reducing Gloves	Size 6 (Small)	0.30mm	5
RR - 2/6.5	Radiation Reducing Gloves	Size 6.5 (Small)	0.30mm	5
RR - 2/7	Radiation Reducing Gloves	Size 7 (Medium)	0.30mm	5
RR - 2/7.5	Radiation Reducing Gloves	Size 7.5 (Medium)	0.30mm	5
RR - 2/8	Radiation Reducing Gloves	Size 8 (Large)	0.30mm	5
RR - 2/8.5	Radiation Reducing Gloves	Size 8.5 (Large)	0.30mm	5
RR - 2/9	Radiation Reducing Gloves	Size 9 (X-Large)	0.30mm	5



a STERIS company

Lumley Close, Thirsk Industrial Park, Thirsk, North Yorkshire, England, YO7 3TD
01845 526660 | info@diagmed.co.uk | diagmed.healthcare

