



## Injection Needles

The widest choice of endoscopic needles available in the UK, designed to ensure safe, effective injection with the greatest ease of operation.

MAKING A DIFFERENCE TO HEALTH

 **DIAGMED**  
HEALTHCARE

a STERIS company

# Carr-Locke Injection Needle

The Carr-Locke Injection Needle was designed to operate in the most tortuous conditions and has become the everyday needle of choice for many around the world.

## The Carr-Locke Injection Needle features:

- Patented design ensures consistent needle projection every time
- Stainless steel, Teflon coated, spring sheath virtually eliminates kinking, whilst offering exceptional functionality and needle performance during challenging procedures
- Ultra-smooth distal metal tip aids in tamponade, and is gentle on friable tissue
- Luer-lock, spring-loaded handle allows for projection and retraction, reducing the risk of inadvertent needle sticks and scope damage
- Requires less than 1cc of fluid to prime (most needles require 1.5- 2.0cc)



## Ordering Information


Product Code	Description	Sheath Diameter (mm)	Length (cm)	Needle Projection	Needle Gauge	Box Size
00711811	Carr-Locke Injection Needle	2.5	230	5mm	25	5
00711812	Carr-Locke Injection Needle	2.5	230	5mm	23	5
00711813	Carr-Locke Injection Needle - Pentax Compatible	2.5	230	5mm	25	5
00711814	Carr-Locke Injection Needle	2.5	230	4mm	23	5

Recommended for EMR and ERCP

“World-renowned for reliability and consistency”

ULTRA SMOOTH METAL HUB  
Aids in tamponade

TEFLON COATED SHEATH  
Allows smooth insertion down endoscope channel



## Articulator™ Injection Needle

The Articulator™ Injection Needle is a reliable needle, ideal for everyday use and dependable for your most immediate needs.

### Features:

From Polypectomy to Oesophageal varices, from ulcers to tattooing, the Articulator™ Injection Needle is the steadfast choice for injection therapy.

- Luer-lock, spring loaded handle to help ensure the needle is engaged or retracted – reducing the risk of inadvertent scope damage and incidental needle sticks, whilst providing maximum needle projection
- Kink resistant, stainless steel sheath to pass smoothly and perform in fully articulated conditions
- Stainless steel hub to aid in tamponade and control needle projection

**ULTRA SMOOTH  
METAL HUB**  
Prevents needle perforations

**STAINLESS STEEL  
SHEATH**  
Kink resistant, smooth passage

“Ideal for everyday use and complicated procedures equally”



### Ordering Information

Product Code	Description	Sheath Diameter (mm)	Length (cm)	Needle Projection (mm)	Needle Gauge	Box Size
00711803	Articulator™ Injection Needle	2.5	160	4	25	5
00711804	Articulator™ Injection Needle	2.5	160	5	25	5
00711807	Articulator™ Injection Needle	2.5	230	5	25	5
00711808	Articulator™ Injection Needle	2.5	230	5	25	5

## Transbronchial Aspiration Needle

The Transbronchial Aspiration Needle is used during a Bronchoscopy to take tissue samples from lymph nodes

## SURECLICK Needle

SURECLICK is an everyday economic needle manufactured in Europe. The 'SURECLICK' needle makes an audible 'click' sound when needle is projected or retracted reducing the risk of inadvertent scope damage and accidental needle injuries.

The process of taking a lymph sample requires the user to undertake a 'blind' procedure where a puncture is made with the needle in the 'narrowing' bronchial canal then pushed into the lymph node.

A syringe (usually 20ml) is attached and then aspiration is performed to draw tissue samples into the needle.


**“Recommended for Bronchoscopy procedures to aspiral fluid/cells to diagnose lung disease”**

### Features:

- Distal metal protector to avoid needle retraction issues
- Teflon coated sheath allows for smooth insertion down endoscope channel
- 23GA is suitable for injection of saline, SPOT, Botox, contrast, adrenaline and other fluids
- 19GA is suitable for injection of Sigmavisc, gelatin based lifting solutions, glue and other thicker solutions
- Economic option

### Ordering Information

Product Code	Description	Sheath Diameter (mm)	Length (cm)	Needle Projection (mm)	Needle Gauge	Box Size
ASP-18-15	Transbronchial Aspiration Needle	1.8	120	15	25	5
09-10-72-32-12AS	SURECLICK	2.3	230	5	23	10
09-11-02-32-12AS	SURECLICK	2.3	230	5	19	10



## Hystoacryl Injection Needle

Application: Needle for the injection of hystoacryl and other viscous solutions.

### Features:

- Robust metal handle for secure push/pull deployment
- Luer lock syringe for secure attachment
- Distal metal tip for greater patient, scope and user safety
- Wide bore needle for greater ease of injection

Recommended  
for variceal  
bleeds



“Ideal for injection of viscous haemostatic agents”

DISTAL METAL TIP

99-09-30-72-32-07  
Hystoacryl Injection  
Needle

### Ordering Information

Product Code	Description	Sheath Diameter (mm)	Length (cm)	Needle Projection (mm)	Needle Gauge	Box Size
99-09-30-72-32-07	Hystoacryl Injection Needle	2.3	230	5	22	5

## EcoClick Injection Needle

Our EcoClick Injection needle is suitable for various Endoscopic injection requirements and offers excellent value for money when compared to other leading brands.

### Features:

- Standard 23-gauge needle size suitable for most Endoscopic use
- 4mm needle projection
- Recommended for SPOT marking and saline injection

### Ordering Information

Product Code	Description	Sheath Diameter (mm)	Sheath Length (cm)	Needle Projection (mm)	Box Size
AG-5082-2423-2304	23ga Series 2 Needle	2.4	230	4	5

“For everyday use with safety in mind”

## MediJet Sclerotherapy Needle

Injection needle for sclerotherapy in the gastrointestinal tract by means of mucosal or submucosal injection.

### Features:

- Locking mechanism which significantly reduces perforation risk
- Tapered distal catheter tip which facilitates atraumatic insertion
- Suitable for gastroscopy and colonoscopy

### Ordering Information

Product Code	Description	Sheath Diameter (mm)	Sheath Length (cm)	Needle Projection (mm)	Box Size
E710230	Medi-Jet 20ga Sclerotherapy Needle	2.3	230	5	10



E710230  
Medi-jet 20ga  
Sclerotherapy needle



AG-5082-2423-2304  
EcoClick Injection  
Needle





Lumley Close, Thirsk Industrial Park, Thirsk, North Yorkshire, England, YO7 3TD  
01845 526660 | [info@diagmed.co.uk](mailto:info@diagmed.co.uk) | [diagmed.healthcare](http://diagmed.healthcare)

