


# Small Channel Devices & Accessories



**DIAGMED**  
HEALTHCARE

a STERIS company



**STERIS**

Endoscopy

# GI Procedure Products

STERIS Endoscopy offers a variety of devices suited for cases requiring endoscopes with smaller channels ranging from 1.0mm to 2.0mm. From foreign body retrieval to tissue sampling and endoscope accessories – we have you covered!


Our specialty bite blocks feature a sturdy, yet forgiving material and smooth insertion rim that help protect delicate instruments and provide comfort for patients.



## Pediatric bite block

**BX00712808**

- » latex-free strap
- » without dental retention rim
- » 38FR



## Intermediate bite block

**BX00712809**

- » latex-free strap
- » without dental retention rim
- » 48FR

## Endo-Boot™

### **BX00711621**

- » single-use foam bumper cushions your endoscope and its delicate components during transport and storage helping to reduce the chance of damage

## ProArmor™ tip protectors

### **BX00711623 / BX00711624**

- » single-use tip protector featuring a flexible, non-porous material and an open design that allows for air flow, reducing the chances of retaining moisture and endoscope damage

## Small Wonder® cleaning brush

### **BX00711616 / BX00711617**

- » specialized, miniature cleaning brush for the smallest endoscope channels (1.00mm - 2.25mm)
- » non-slip, braided wire shaft allows for steady grip


# Polypectomy and Tissue Acquisition

World renowned for its reliability and consistency, the Carr-Locke injection needle is a device that can be trusted during everyday procedures as well as difficult maneuvers.

## Carr-Locke injection needle – mini

### **BX00711822**

- » designed to operate in tortuous conditions, including duodenoscopes, and is ideal for injection therapy across a broad range of smaller channel endoscopes
- » patented design ensures consistent needle projection every time



STERIS Endoscopy's line of forceps devices address a variety of needs for tissue sampling. Our disposable biopsy forceps feature a Teflon sheath for ease of insertion and are compatible with pediatric bronchoscopes and gastroscopes.

## Disposable biopsy forceps

### **BX00711261 / BX00711262**

- » produced with sharpened cutting edges
- » built with stainless steel cups and state-of-art technology
- » disposable forceps are available with or without a spike



# Foreign Body Management

Foreign body retrieval can be challenging. Our unique line of products specifically addresses these difficult procedures.



## Raptor® grasping device - mini

### **BX00711178**

- » combines alligator and rat tooth capabilities into one device
- » smaller hybrid jaw and catheter ideal for retrieval of foreign bodies in tight positions or procedures requiring an ultra-slim endoscope

## Roth Net® retriever - mini

### **BX00711057 / BX00711168**

- » specifically designed for retrieval of a variety of foreign bodies or polyps in smaller channel endoscopes
- » unique design maintains shape and opens to full capacity across multiple passes



# Product Information

Procedural Infection Prevention						
part number	description	sheath diameter (mm)	channel size (mm)	bristle diameter(mm)	length (cm)	qty/box
BX00711616	Small Wonder cleaning brush	0.5	1.00-2.25	2.25	160	15
BX00711617	Small Wonder cleaning brush	0.5	1.00-2.25	2.25	230	15
part number	description		compatible scope OD (mm)			qty/box
BX00711621	small endoscope tip protector		3.0 - 8.8			50
BX00711623	ProArmor endoscope tip protector - sterile		4.9-15			50
BX00711624	ProArmor endoscope tip protector - non-sterile		4.9-15			50
part number	description		opening size	dental retention rim		qty/box
BX00712808	pediatric bite block		12.6mm / 38FR	no		25
BX00712809	intermediate bite block		16mm / 48FR	no		100
Polypectomy & Tissue Acquisition						
part number	description		sheath diameter(mm)	length (cm)	qty/box	
BX00711261	disposable biopsy forceps - oval with spike (sterile)		1.8	160	5	
BX00711262	disposable biopsy forceps - oval (sterile)		1.8	160	5	
BX00711822	Carr-Locke injection needle - 5mm projection / 25 gauge (sterile)		1.8	230	5	
Foreign Body Management						
part number	description	sheath diameter (mm)	approx. net size (cm)	length (cm)	qty/box	
BX00711057	Roth Net retriever - mini	1.8	160	2 x 4.5	5	
BX00711168	Roth Net retriever - mini (sterile)	1.8	160	2 x 4.5	5	
BX00711178	Raptor grasping device - mini (sterile)	1.9	200	-	5	



a STERIS company

Lumley Close, Thirsk  
Industrial Park, Thirsk,  
North Yorkshire,  
England, YO7 3TD  
01845 526660 |  
[info@diagmed.co.uk](mailto:info@diagmed.co.uk) |  
[diagmed.healthcare](http://diagmed.healthcare)

© 2020 STERIS Endoscopy or its affiliates. All rights reserved.  
Unless otherwise indicated, all marks denoted with ® or ™ are registered with the  
U.S. Patent and Trademark Office, or trademarks owned by STERIS or its affiliates.

Follow us:    